




# New Parts to fit Hill-Rom Stabilet Infant Warmer

RPI Part #AIB116  
OEM Part #BDZ107/BDZ108

*RPI Added Value!*

## FLUORESCENT BALLAST

• Generic part; can be used to replace right hand or left hand ballast  
**Fits:** Overhead Light Assy  
**Models:** Stabilet



RPI Part #AIH118  
OEM Part #BCR730

## FLUORESCENT TUBE HOLDER


**Fits:** Overhead Light Assy.  
**Models:** Stabilet



RPI Part #AIE115  
OEM Part #BDJ002 & BDJ003

## ENCLOSURE, FLUORESCENT TUBE


**Fits:** Overhead Light Assy.  
**Models:** Stabilet



RPI Part #LMP021  
OEM Part #BDJ001

## FLUORESCENT LAMP

**Fits:** Overhead Light Assy  
**Models:** Stabilet




RPI Part #AIR112  
OEM Part #BDZ205

*RPI Added Value!*

## LIGHT RELAY ASSY


**Fits:** Overhead Light Assy.  
• Includes installed timing resistor  
**Models:** 2000A/B, 2000C, 1250



RPI Part #AIR124  
OEM Part #BCZ952

## RECEPTACLE ASSY, PROBE


**Fits:** Chassis Assy  
**Models:** Stabilet



RPI Part #AIB109  
OEM Part #BKK010

## BASSINET BRACKET/HINGE

• Beige  
**Fits:** Bassinet  
**Models:** Stabilet 2000A/B, 2000C, 1250



RPI Part #AIB110  
OEM Part #KB010

## BASSINET BRACKET/HINGE

• Black  
**Fits:** Bassinet  
**Models:** Stabilet 2000A/B, 2000C, 1250



RPI Part #AIB111  
OEM Part #BKK015

## CORNER BUMPER


**Fits:** Base Assy.  
**Models:** 2000A/B, 2000C



RPI Part #AIC119  
OEM Part #BPM002

## CASTER - 5" SWIVEL WITH BRAKE


**Fits:** Base Assembly  
**Models:** Stabilet



RPI Part #AIC120  
OEM Part #BPM002

## CASTER - 5" SWIVEL (NO BRAKE)


**Fits:** Base Assembly  
**Models:** Stabilet



RPI Part #AIS113  
OEM Part #BC0003

## THERMAL SWITCH


**Fits:** Overhead Light Assy.  
**Models:** 2000A/B, 2000C, 1250



RPI Part #AIT114  
OEM Part #BCP425 & BEZ504

## FAST PAK TRIAC ASSY WITH CABLE


**Fits:** Chassis Assy.  
**Models:** 2000A/B, 2000C, 1250



RPI Part #RPC464  
OEM Part # NA

## HEAT SINK COMPOUND

• MSDS and Technical Data sheet available on RPI website  
• Non-Silicone based thermal grease  
• 5 per package  
• Size: .07 fl. oz (2ml)




RPI Part #AIS122  
OEM Part #BDB001

*RPI Added Value!*

## SWITCH, MAIN POWER


• Includes AIL123 and LMP022  
**Fits:** Front Panel  
**Models:** Stabilet 2000C



RPI Part #AIL123  
OEM Part #BDB002

## LENS, MAIN POWER SWITCH


**Fits:** Main Power Switch - RPI Part #AIS122  
**Models:** Stabilet 2000C



RPI Part #LMP022  
OEM Part #BDB003

## LAMP

• 14 VDC, 1.1W, Base: Wedge, T-1 3/4  
• Incandescent  
**Fits:** Main Power Switch - RPI Part #AIS122  
**Models:** Stabilet 2000C



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Hill-Rom Air Shields equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 02/05