




New parts to fit Pelton & Crane Magnaclave Sterilizers!

NEW




RPI Part #PCR251
OEM Part #045805

RELAY (K1 & K3)

- Size: 2.69" L x 1.32" W x 1.85" H
- Coil: 200/240 VAC
- Contacts: DPST; 25A N.O.
- 1/4" Tabs

Fits: Chassis
Model: Magnaclave
(Serial #7525 and above)

NEW




RPI Part #PCR252
OEM Part #045806

RELAY (K2)

- Size: 2.70" L x 1.36" W x 1.49" H
- Coil: 220/240 VAC
- Contacts: DPDT; 30A N.O./3A N.C.
- 1/4" Tabs

Fits: Chassis
Model: Magnaclave
(Serial # 7525 and above)

NEW




RPI Part #PCR234
OEM Part #041002 & A62-011081

RELAY (ELECT. CONTACTOR)

- Coil: 208/240VAC
- 2 pole, contacts rated 20A
- 1/4" Tabs

Fits: Rear Chassis Support
Model: Magnaclave

NEW




RPI Part #PCK245
OEM Part #004536 (Knob only) & 004507 (Faceplate only)

KNOB (TIMER)

- Knob: Black
- Time: 0-60 minutes
- Fits 1/4" dia. shaft
- Includes set screw
- Also available: Electric Timer (60 Minute) (RPI Part #PCT081)

Fits: Electric Timer (60 Minute)
Model: Magnaclave

NEW




RPI Part #PCT081
OEM Part #004531 & 041001

ELECTRIC TIMER (60 MINUTE)

- DPDT; 60 minute
- Switch: 7.5A @ 240VAC; 15A @ 125VAC
- Includes mounting hardware
- Also available: Buzzer (RPI Part #PCB088) and Knob (Timer) (RPI Part #PCK245)

Fits: Control Panel
Model: Magnaclave

NEW




RPI Part #PCB088
OEM Part #004491

BUZZER

- 5mA @ 220VAC, 60 Hz.
- 1/4" tabs
- Also available: Electric Timer (60 Minute) (RPI Part #PCT081)

Fits: Mounts on machine
Model: Magnaclave

NEW




RPI Part #PCK243
OEM Part #004625 (Knob only), 004508 (Faceplate only) & 004626 (Stop Plate only)

KNOB (TEMPERATURE)

- Knob: Black
- Faceplate: 240°F - 270°F
- Fits 1/4" dia. shaft
- Includes set screw

Fits: Thermostat Control
Model: Magnaclave

NEW



RPI Part #PCK244
OEM Part #004536 (Knob only) & 004506 (Faceplate only)

KNOB (FUNCTION)

- Knob: Black
- Faceplate: Power Off, Fill, Sterilize and Vent
- Fits 1/4" dia. shaft
- Includes set screw

Fits: Main Valve Control Shaft
Model: Magnaclave

Also available,
new parts to fit
Pelton & Crane
OCM, OCR & OCR+
Table Top Sterilizers
-see back for details



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 02/12



replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD
U.S. POSTAGE PAID
VAN NUYS, CA
PERMIT NO. 2399

ADDRESS SERVICES REQUESTED


**New parts to fit
Pelton & Crane
table top sterilizers!**

Model: OCM, OCR & OCR+



New parts to fit Pelton & Crane OCM, OCR & OCR+ Sterilizers!

RPI Part #PCK246
OEM Part #010784 (Complete Knob),
013016 (Knob only) & **004178** (Skirt only)




KNOB (THERMOSTAT)

- Knob: Black
- Skirt: Stainless steel
- Fits: 1/4" dia. shaft
- Includes set screw
- Also available: Controller (RPI Part #PCC112)

Fits: Controller Potentiometer

Models: OCM , OCR and OCR+

RPI Part #PCK248
OEM Part #013015



KNOB (TIMER & FUNCTION)

- Knob: Black
- Fits 1/4" dia. shaft
- Includes set screw
- Also available: Bell Timer (RPI Part #PCT049)

Fits: Timer and Main Valve Control Shaft

Models: OCM, OCR, OCR+


RPI Part #PCC112
OEM Part #004144

CONTROLLER

- Includes mounting hardware
- Used with Solid State Controller
- 20 K Ω
- .250" dia. shaft
- Also Included in: Thermostat Conversion Kit (RPI Part #PCK128)
- Also available: Knob (Thermostat) (RPI Part #PCK246),
Solid State Controller with Harness (RPI Part #PCT007) and
Solid State Controller (RPI Part #PCT169)

Fits: Control Panel

Models: OCM, OCR, OCR+ & Sentry




RPI Part #PCT049
OEM Part #004112

BELL TIMER

- Time & Alert Specs: 60 minutes, Bell Alert
- Mount Specs: 3/8-24 x 9/32" lg. Bushing
- Shaft: 1/4" dia. x 1/2" lg.
- Size: 2-11/16" dia. x 1-1/4" lg.
- Includes mounting hardware

Fits: Control Panel

Models: OCM, OCR, OCR+ and Sentry



**Also available, new parts to fit Pelton & Crane
Magnaclave Table Top Sterilizers - see other side for details**



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 02/12