



Table Top Sterilizers

NEW PARTS TO FIT TUTTNAUER

FITS 1730, 2340, 2540, 3140, 3850 & 3870 SERIES


RPI Part #TUK075

OEM Part # (No OEM # Available)

RESERVOIR COVER WITH DIPSTICK

- Includes all parts as shown

Fits: Water Reservoir/Top Outer Cover



RPI Part #TUB072

OEM Part # (No OEM # Available)

RUBBER BOOT FOR DOOR SWITCH

- This push button rubber boot will extend the life of the door switch by preventing moisture from entering into it

Fits: Door Switch




RPI Part #TUS014

OEM Part #01910190

DOOR SWITCH

- 15A @ 125-250V
- Includes mounting nut

Fits: Front panel



RPI Part #TUB064

OEM Part #CT245010

DOOR BELLOWS HOUSING BOLT

- 3/8-19 BSPP; Brass

Fits: Door




RPI Part #TUT073

OEM Part # TEST-3

TAP (3/8-19 BSPP)

- Use to re-tap bellows housing in door



RPI Part #TUG022

OEM Part #02610029

WATER RESERVOIR GASKET

- ~6" round bellows style
- Silicone; Color: black

Fits: Water Reservoir



RPI Part #TUC067

OEM Part #02550019

WATER RESERVOIR/FILTER COVER

- Cover for both the Water Reservoir **and** the HEPA Filter compartment

Fits: Water Reservoir **and** HEPA Filter compartment



RPI Part #TUS068

OEM Part #02550043

WATER RESERVOIR DIPSTICK

Fits: Water Reservoir Cover




RPI Part #TUV065

OEM Part #03110002

SAFETY VALVE (40 PSI)

- Brass; 1/4" MPT
- Cracking pressure: 40 PSI

Fits: Safety Valve Holder



40 PSI




RPI Part #TUV011

OEM Part #03110003

SAFETY VALVE (37 PSI)

- Brass; 1/4" MPT
- Cracking pressure: 37 PSI

Fits: Safety Valve Holder



37 PSI



RPI Part #MIF062


OEM Part # (No OEM # Available)

FILTERS (mesh style)

- 2 per package
- Size: 30 x 30 mesh (stainless steel)

- Can be used to filter chamber debris at fill and vent ports

Fits: Chamber ports




RPI Part #RPR583

OEM Part #02819993

AC INLET RECEPTACLE

- 15A@250V
- Contact Style: .250" x .032" (6.3mm x 0.8mm) tab





RPI Part #RPC582

OEM Part #02819996

POWER CORD

- 14/3 SJT, 105°C; 15A@220V
- Plug Type: NEMA 6-15P
- Connector Type: EN60320-C13R Right Angle

Fits: AC Inlet Receptacle





RPI Part #TUC028

OEM Part #02819995

POWER CORD

- 14/3 SJT, 105°C; 15A@125V
- Plug Type: NEMA 5-15P
- Connector Type: EN60320-C13R Right Angle

Fits: AC Inlet Receptacle



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com


The above parts are manufactured by Replacement Parts Industries, Inc. to fit Tuttnauer equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer WEB 05/07

CONDENSATION COILS ... RPI HAS THE ONE YOU NEED!

Over the years, Tuttnauer has changed the configuration of their reservoirs by reversing the coil.
RPI has you covered – we offer all three styles!


RPI Part #TUC040
OEM Part #CU836101

Fits 1730E/EK
& 1730M/MK




RPI Part #TUC041
OEM Part #CT836101

Fits original 2340 & 2540 Series
Coil joins reservoir at the
left rear and vents
towards the front
of the machine



RPI Part #TUC063
OEM Part #CT836101

Fits newer 2340 & 2540 Series
Coil joins reservoir at the
left front and vents
towards the rear
of the machine



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028
EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

THIS CHART LISTS ALL OF THE RPI PARTS (AS OF MAY, 2007) AVAILABLE TO FIT TUTTNAUER EQUIPMENT

RPI Part #	OEM Part #	Description	Model 1730E	Model 1730EK	Model 1730M	Model 1730MK	Model 2340E	Model 2340EK	Model 2340M	Model 2340MK	Model 2540E	Model 2540EK	Model 2540M	Model 2540MK	Model 3870E	Model 3870EK	Model 3870M	Model 3870MK
MDL019	01960007	SIGNAL LIGHT - RED																
MDL021	01910257	SIGNAL LIGHT - ORANGE																
MDL022	01910258	SIGNAL LIGHT - GREEN																
MIF062	(No OEM Part # Available)	FILTERS (mesh style)																
RPA459	(No OEM Part # Available)	PIPE SEALANT 567																
RPC582	02819996	POWER CORD																
RPH285	02620016A	DRAIN HOSE																
RPO386	(No OEM Part # Available)	O-RINGS																
RPO387	(No OEM Part # Available)	O-RING																
RPR583	02819993	AC INLET RECEPTACLE																
TUB019	(No OEM Part # Available)	DOOR BELLOWS																
TUB023	01910098	CIRCUIT BREAKER		1		1	1	2	1	2	1	2	1	2	2		2	
TUB024	01910097	CIRCUIT BREAKER	3	2	3	2	2		2	2		2						
TUB047	01910100	CIRCUIT BREAKER (15A)	1	3	1	3	3	2	3	2	3	2	3	2		2		
TUB048	01910099	CIRCUIT BREAKER (7A)	2		2													
TUB064	CT245010	DOOR BELLOWS HOUSING BOLT																
TUB072	(No OEM Part # Available)	RUBBER BOOT FOR DOOR SWITCH																
TUC028	02819995	POWER CORD																
TUC040	CU836101	CONDENSATION COIL																
TUC041	CT836101	CONDENSATION COIL																
TUC063	CT836101	CONDENSATION COIL																
TUC067	02550019	WATER RESERVOIR/FILTER COVER																
TUG001	02610020	DOOR GASKET		1														
TUG002	02610005	DOOR GASKET						6										
TUG003	02610023	DOOR GASKET																
TUG012	02300011	PRESSURE GAUGE WITH																
TUG021	02610118	DOOR GASKET						7	7									
TUG022	02610029	WATER RESERVOIR GASKET																
TUH004	01720011	HEATER ELEMENT			1													
TUH005	01720002	HEATER ELEMENT																
TUH006	01720003	HEATER ELEMENT																
TUH015	(No OEM Part # Available)	ATTACHING HARDWARE																
TUH016	01720012	HEATER ELEMENT		2		2												
TUH017	01720013	HEATER ELEMENT																
TUH018	01720014	HEATER ELEMENT																
TUH027	01720001	HEATER ELEMENT																
TUH031	CU841010	SAFETY VALVE HOLDER																
TUH032	CT841020	SAFETY VALVE HOLDER																
TUH043	CU900012	WIRE HARNESS			5	5												
TUH044	CT900012	WIRE HARNESS						5	5			5	5					
TUJ033	CB842010	AIR VALVE JET																
TUJ034	CT842010	AIR VALVE JET																
TUK030	CT241111	DOOR BELLOWS ASSEMBLY KIT																
TUK037	(No OEM Part # Available)	MULTI-PURPOSE VALVE REPAIR KIT																
TUK049	02450002	TIMER & MULTI-PURPOSE VALVE KNOB																
TUK050	02450003	THERMOSTAT KNOB																
TUK053	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT																
TUK054	(No OEM Part # Available)	SAFETY VALVE HOLDER KIT																
TUK055	N/A	MOUNTING HARDWARE																
TUK075	(No OEM Part # Available)	RESERVOIR COVER WITH DIPSTICK																
TUS009	01910172	POWER SWITCH																
TUS013	01910191	MICROSWITCH																
TUS014	01910190	DOOR SWITCH																
TUS035	01620301	OVER TEMPERATURE SAFETY SWITCH		4		4		4		4		4		4				
TUS057	01910197	MICRO-SWITCH (NEW STYLE)																
TUS068	02550043	WATER RESERVOIR DIPSTICK																
TUT007	01620101	CONTROL THERMOSTAT																
TUT008	01910011	TIMER (Old OEM Part #01910005)																
TUT038	01620004	SAFETY THERMOSTAT (MANUAL RESET)																
TUT039	01620103	SAFETY THERMOSTAT (AUTOMATIC RESET)																
TUT073	TEST-3	TAP (3/8-19 BSPP)																
TUV011	03110003	SAFETY VALVE (37 PSI)																
TUV025	CT810013	MULTI-PURPOSE VALVE ASSEMBLY																
TUV042	CT844180	DRAIN VALVE ASSEMBLY																
TUV065	03110002	SAFETY VALVE (40 PSI)																

PARTS AVAILABLE TO FIT TUTTNAUER EQUIPMENT

1 Fits 110V units 2 Fits 220V units 3 Fits 110V and 220V units 4 Only models prior to 1993 5 Existing connections should accommodate most units manufactured after 01/93. However some modifications may be necessary.
6 Fits units with S/N 8805 and below 7 Fits units with S/N 8806 and above