



# NEW! PCB's & SWITCHES TO FIT PELTON & CRANE CHAIRMAN, CHAIRMAN 5000 & COACHMAN DENTAL CHAIRS



## CHAIRMAN


RPI Part #PCS720  
New OEM Part #007899

### EXIT SWITCH ASSEMBLY

- Normally Open
- 3/4A @ 125VAC
- Black button

**Fits:** Chair back assembly

**Models:** Chairman




RPI Part #PCS719  
New OEM Part #007898

### OPERATE SWITCH ASSEMBLY

- Normally Open
- 3/4A @ 125VAC
- Red button

**Fits:** Chair back assembly

**Models:** Chairman




RPI Part #PCS713  
OEM Part #007431

### ROLLER MICRO SWITCH

- 5A @ 250 VAC

**Fits:** Seat wiring assembly

**Model:** Chairman




RPI Part #PCS714  
New OEM Part #007506

### FORMED LEVER SWITCH

- 5A @ 120VAC

**Fits:** Base, lift & attachment assembly

**Models:** Chairman




RPI Part #PCS715  
OEM Part #007596

### TRAVERSE LIMIT SWITCH

- 5A @ 125/250VAC

**Fits:** Base, lift & attachment assembly

**Models:** Chairman




RPI Part #PCP721  
OEM Part #007595

### FOOT SWITCH PEDAL

- Black plastic

**Fits:** Chair base foot switch assembly

**Models:** Chairman




RPI Part #PCS716  
OEM Part #007429

### TRAVERSE SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (dark grey paddle)

**Fits:** Back assembly and foot switch assembly

**Model:** Chairman




RPI Part #PCS717  
OEM Part #007428

### TILT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V; Snap in mounting with low profile paddle (white paddle)

**Fits:** Back assembly (Chairman) and Foot Switch Assembly (Chairman 5090)

**Model:** Chairman & Chairman 5090




RPI Part #PCS718  
OEM Part #007427

### LIFT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (black paddle)

**Fits:** Back assembly

**Model:** Chairman



## COACHMAN & COACHMAN II

RPI Part #PCC728  
OEM Part #94-416-43

### LIFT CAPACITOR (56-60 µf)

- Plastic; Electrolytic medium

**Fits:** Base motor/lift motor

**Models:** Chairman 5000 & Coachman



RPI Part #PCC734  
OEM Part #77-31-904/019518

### CAPACITOR (21-25 µf)

- Plastic; Electrolytic medium

**Fits:** Incline/recline and upper traverse

**Models:** Chairman 5000 & Coachman




RPI Part #PCB723  
New OEM Part #018326

### SWITCH BACK PCB ASSEMBLY

**Fits:** Chair back assembly

**Models:** Coachman




RPI Part #PCS724  
OEM Part #018420

### MICRO SWITCH

- 5A @ 250 VAC

**Fits:** Chair back assembly

**Model:** Coachman




RPI Part #PCS722  
OEM Part #023944

### FOOT SWITCH MICROSWITCH

- 15A @ 125VAC

**Fits:** 4-way foot switch assembly

**Models:** Coachman & Coachman II




RPI Part #PCS725  
OEM Part #018956

### MOTOR LIMIT SWITCH

- 15A @ 120VAC

**Fits:** Base assembly

**Models:** Coachman & Coachman II



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 11/04



# NEW! PCB's & SWITCHES TO FIT PELTON & CRANE CHAIRMAN, CHAIRMAN 5000 & COACHMAN DENTAL CHAIRS



## CHAIRMAN 5000


**RPI Part #PCC734**  
**OEM Part #77-31-904/019518**

**CAPACITOR (21-25 µf)**  
• Plastic, Electrolytic medium  
**Fits:** Incline/recline and upper traverse  
**Models:** Chairman 5000 and Coachman




**RPI Part #PCC728**  
**OEM Part #94-416-43**

**LIFT CAPACITOR (56-60 µf)**  
• Plastic; Electrolytic medium  
**Fits:** Base motor/lift motor  
**Models:** Chairman 5000 & Coachman




**RPI Part #PCB736**  
**New OEM Part #94 34 424**

**UP INTERCONNECT (PCB)**  
**Fits:** Incline/recline controls  
**Models:** Chairman 5000 & 5010




**RPI Part #PCS735**  
**OEM Part #77 40 947**

**A/C POWER SWITCH**  
• 16A @ 125/250VAC  
• Use only with 240V models  
**Fits:** Inner base assembly  
**Models:** Chairman 5000 (**240V models only**)



**RPI Part #PCS731**  
**OEM Part #15-25-737**


**CORRECTION LIMIT SWITCH**  
• 10A @ 125 VAC  
**Fits:** Lift motor  
**Model:** Chairman 5000



## CHAIRMAN 5010


**RPI Part #PCB736**  
**New OEM Part #94 34 424**

**UP INTERCONNECT (PCB)**  
**Fits:** Incline/recline controls  
**Models:** Chairman 5000 & 5010



**RPI Part #RPF367**  
**OEM Part #18-81-163**


**FUSE (2A, 250V)**  
• 5 per package  
• Size: 5 mm x 20 mm (.197" x .796")  
• MDL Glass, Fast-Acting  
**Fits:** Inner base assembly (**240V models only**)  
**Models:** Chairman 5000 & 5010



## CHAIRMAN 5090

**RPI Part #PCB733**  
**New OEM Part #30 06 629**

**LIMIT SWITCH (PCB)**  
**Fits:** Incline/recline and upper traverse  
**Models:** Chairman 5090




**RPI Part #PCA737**  
**OEM Part #18 99 702**

**FOOT SWITCH & CABLE ASSEMBLY**  
• Momentary contact (On/Off/On)  
• 6A @ 125 V; 3A @ 250 V; Snap in mounting  
with low profile paddle (white paddle)  
**Fits:** Inner base  
**Models:** Chairman 5090



**RPI Part #PCS717**  
**OEM Part #007428**

**TILT SWITCH**  
• Momentary contact (On/Off/On)  
• 6A @ 125 V; 3A @ 250 V; Snap in mounting  
with low profile paddle (white paddle)  
**Fits:** Back assembly (Chairman) and Foot Switch  
Assembly (Chairman 5090)  
**Model:** Chairman & Chairman 5090



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 11/04