





New Parts to Fit

AIR TECHNIQUES DENTAL FILM PROCESSORS


A/T 2000®, A/T 2000® Plus & A/T 2000® XR




RPI Part #ATA670 (Plastic mounting plate not included)
OEM Part #43330
HEATER AND FAN ASSEMBLY
Fits: Cover Assembly
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR




RPI Part #ATH668
OEM Part #43811
HEATER ELEMENT
Fits: Heater and Fan Assembly
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR




RPI Part #ATF667
OEM Part #43530
FAN
Fits: Heater and Fan Assembly
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR



RPI Part #ATT669
OEM Part # (No OEM # Available)
OVERHEAT THERMOSTAT
Fits: Heater Element
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR



RPI Part #RPH638
OEM Part #90219
FUSE HOLDER (MAIN)
Fits: Back panel
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR



RPI Part #RPF552
OEM Part #41914
FUSE (8A, 250V) - TIME DELAY - 5/PKG
Fits: Fuse Holder
Models: A/T 2000, A/T 2000 Plus & A/T 2000 XR

Listing of RPI Parts to fit Air Techniques A/T 2000, A/T 2000 Plus & A/T 2000 XR Dental Film Processors

RPI PART #	OEM PART #	DESCRIPTION	A/T 2000	A/T 2000 Plus	A/T 2000 XR
ATA644	43224-1	REMOTE SOLENOID WATER VALVE ASSEMBLY	•	•	
ATA670	43330	HEATER AND FAN ASSEMBLY	•	•	•
ATB643	43420	BASE PC BOARD	•		
ATF667	43530	FAN	•	•	•
ATH639	43290	HEATER BAR ASSEMBLY	•	•	
ATH647	43440	BASE WIRING HARNESS	•		
ATH668	43811	HEATER ELEMENT	•	•	•
ATK646	43942	DRIVE MOTOR BOARD WITH HARNESS	•		
ATS641	43447	WATER LEVEL SENSOR	•	•	
ATS642	30552	THUMB SCREW - 8/PKG	•	•	
ATT669	(No OEM # Available)	OVERHEAT THERMOSTAT	•	•	•
ATV645	(No OEM # Available)	SOLENOID VALVE ASSEMBLY	•	•	
ATW638	43271	DRAIN TUBE WASHER - 6/PKG	•	•	•
RPF552	41914	FUSE (8A, 250V) - TIME DELAY - 5/PKG	•	•	•
RPH638	90219	FUSE HOLDER (MAIN)	•	•	•
RP0244	43362	O-RING - 6/PKG	•	•	



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Air Techniques equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Web Flyer 09/08

DENTAL FILM PROCESSORS – AIR TECHNIQUES

SECTION 4

PAGE 601A

UPDATED 09/08