











ADAPTER AND LEADWIRE ACCESSORIES

RPI PART #	DESCRIPTION	ACCESSORY	DETAILS
LAC002	3 LEAD COMBINER		<ul style="list-style-type: none"> • Number of Leads: 3 • Use with Dual/Shielded Leadwires only • Included in: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
LAC003	5 LEAD COMBINER		<ul style="list-style-type: none"> • Number of Leads: 5 • Use with Dual/Shielded Leadwires only • Included in: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)
LAK001	LEADWIRE ADAPTER KIT	  <p>Includes 10 adapters individually labeled with leadwire identifiers</p>	<ul style="list-style-type: none"> • Includes 10 adapters individually labeled with leadwire identifiers • Adapts banana plugs or pins to tab or snap electrodes • Widest contact area for superior trace quality • Flat, wide bodied design helps prevent twisting • Durable construction - built to last
LAK006	3 LEAD SHIELDED ACCESSORY KIT	 	<ul style="list-style-type: none"> • Number of Leads: 3 • Use with Dual/Shielded Leadwires only • Includes: <ul style="list-style-type: none"> 3 Lead Combiner (RPI Part #LAC002) and 3 Lead Organizer (RPI Part #LAR004)
LAK007	5 LEAD SHIELDED ACCESSORY KIT	 	<ul style="list-style-type: none"> • Number of Leads: 5 • Use with Dual/Shielded Leadwires only • Includes: <ul style="list-style-type: none"> 5 Lead Combiner (RPI Part #LAC003) and 5 Lead Organizer (RPI Part #LAR005)
LAR004	3 LEAD ORGANIZER		<ul style="list-style-type: none"> • Number of Leads: 3 • Included in: 3 Lead Shielded Accessory Kit (RPI Part #LAK006)
LAR005	5 LEAD ORGANIZER		<ul style="list-style-type: none"> • Number of Leads: 5 • Included in: 5 Lead Shielded Accessory Kit (RPI Part #LAK007)