

















REUSEABLE: DIN/UNSHIELDED LEADWIRE SETS

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LKM002	3 LEADWIRE SET - 24" DIN/PINCH	3	~24"	Din/Unshielded .060" Safety Socket 	Pinch  WHITE RED BLACK	Also available: • Disposable Leadwire Set (RPI Part #LDM308) & 3 Lead Organizer (RPI Part #LAR004)
LKM005	3 LEADWIRE SET - 24" DIN/SNAP	3	~24"	Din/Unshielded .060" Safety Socket 	Snap  WHITE RED BLACK	Also available: • Disposable Leadwire Set (RPI Part #LDM301) & 3 Lead Organizer (RPI Part #LAR004)
LKM008	3 LEADWIRE SET - 24" DIN/PINCH	3	~24"	Din/Unshielded .060" Safety Socket 	Pinch  WHITE GREEN BLACK	Also available: • Disposable Leadwire Set (RPI Part #LDM309) & 3 Lead Organizer (RPI Part #LAR004)
LKM011	3 LEADWIRE SET - 24" DIN/SNAP	3	~24"	Din/Unshielded .060" Safety Socket 	Snap  WHITE GREEN BLACK	• This Leadwire Set replaces RPI Part #HPC003 Also available: • Disposable Leadwire Set (RPI Part #LDM302) & 3 Lead Organizer (RPI Part #LAR004)
LKM026	5 LEADWIRE SET - 24" DIN/PINCH	5	~24"	Din/Unshielded .060" Safety Socket 	Pinch  RED WHITE GREEN BLACK BROWN	Also available: • Disposable Leadwire Set (RPI Part #LDM310) & 5 Lead Organizer (RPI Part #LAR005)
LKM029	5 LEADWIRE SET - 24" DIN/SNAP	5	~24"	Din/Unshielded .060" Safety Socket 	Snap  RED WHITE GREEN BLACK BROWN	Also available: • Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)
LKM030	5 LEADWIRE SET - 40" DIN/SNAP	5	~40"	Din/Unshielded .060" Safety Socket 	Snap  RED WHITE GREEN BLACK BROWN	• This Leadwire Set replaces RPI Part #HPC007 Also available: • 24" Disposable Leadwire Set (RPI Part #LDM303) & 5 Lead Organizer (RPI Part #LAR005)

REUSEABLE: DIN/UNSHIELDED LEADWIRE SETS *(Continued)*

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	NOTES
LKM038	7 LEADWIRE SET - 24" DIN/PINCH	7	~24"	Din/Unshielded .060" Safety Socket	Pinch 	Also available: <ul style="list-style-type: none"> Disposable Leadwire Set (RPI Part #LDM311), 3 Lead Organizer (RPI Part #LAR004), 5 Lead Organizer (RPI Part #LAR005)
LKM041	7 LEADWIRE SET - 24" DIN/SNAP	7	~24"	Din/Unshielded .060" Safety Socket	Snap 	Also available: <ul style="list-style-type: none"> Disposable Leadwire Set (RPI Part #LDM304), 3 Lead Organizer (RPI Part #LAR004), 5 Lead Organizer (RPI Part #LAR005)