









# New Parts to fit Castle/Getinge® Bulk Sterilizers

<p><b>Gasket shapes to fit door!</b></p> <p><b>NEW</b></p> <p><b>RPI Part #CSG030</b> <b>OEM Part #513966</b></p> <p><b>DOOR GASKET (16" x 16")</b></p> <ul style="list-style-type: none"> <li>• Material: blue silicone</li> <li>• 20.50" OD</li> <li>• "U" shaped cross section</li> </ul> <p><b>Fits: 16" x 16" Door</b> <b>Models: 3522, 3523, 3525, 3533, 3622, 3623 &amp; 3633</b></p>   <p><small>"U" SHAPED CROSS SECTION SHOWN ACTUAL SIZE</small></p>	<p><b>Gasket shapes to fit door!</b></p> <p><b>NEW</b></p> <p><b>RPI Part #CSG032</b> <b>OEM Part #536651</b></p> <p><b>DOOR GASKET (24" x 36")</b></p> <ul style="list-style-type: none"> <li>• Material: orange silicone</li> <li>• 38.00" OD</li> <li>• 1" OD round cross section</li> <li>• The "Big Orange" gasket</li> </ul> <p><b>Fits: 24" x 36" Door</b> <b>Models: MC3600</b></p>   <p><small>1" OD ROUND CROSS SECTION SHOWN ACTUAL SIZE</small></p>
<p><b>Gasket shapes to fit door!</b></p> <p><b>NEW</b></p> <p><b>RPI Part #CSG031</b> <b>OEM Part #513967</b></p> <p><b>DOOR GASKET (20" x 20")</b></p> <ul style="list-style-type: none"> <li>• Material: blue silicone</li> <li>• 25.50" OD</li> <li>• "U" shaped cross section</li> </ul> <p><b>Fits: 20" x 20" Door</b> <b>Models: 3522, 3523, 3525, 3533, 3622, 3623 &amp; 3633</b></p>   <p><small>"U" SHAPED CROSS SECTION SHOWN ACTUAL SIZE</small></p>	

**NEW**

**SOL-1 & 2**  
**Steam to Chamber & Fast Exhaust**


**RPI Part #CSK024**  
**OEM Part #61301604084**

**SOLENOID VALVE REPAIR KIT**

- Diaphragm style repair kit
- EPDM seals

**Fits: 1" Diaphragm Style Valve**  
**Models:**

- 3622, 3623 & 3633 (Sol-1 Steam to Chamber)
- 3622, 3623 & 3633 (Sol-2 Fast Exhaust)









# replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

LIST OF RPI PARTS TO FIT CASTLE/GETINGE®  
AS OF FEBRUARY 2010.

RPI PART #	OEM PART #	DESCRIPTION	MODELS
AMK174	56161	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533
AMK177	61301601739	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533, 3622, 3623, 3633
AMK183	54190	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533
CSC063	(No OEM Part # Available)	Coil	3522, 3533, 3622, 3623, 3633
CSC064	(No OEM Part # Available)	Coil	3523
CSD051	517153	Diaphragm Replacement 1.0	3522, 3523, 3525, 3533
CSF052	61301600303	Filter, Air-In	100, 200, 400/500, 3522, 3523, 3525, 3533, 3622, 3633, 4200MPS, 6900AMPS
CSG001	61301605016 58398 (old)	Door Gasket (20" X 20")	100, 800, 3522, 3523, 3533
CSG002	61301605017 58397 (old)	Door Gasket (16" X 16")	100, 3522, 3523, 3525, 3533
CSG003	40441	Door Gasket (Green)	6300, 6400, 6900
 <b>CSG030</b>	513966	Door Gasket (16" X 16")	3522, 3523, 3525, 3533, 3622, 3623, 3633
 <b>CSG031</b>	513967	Door Gasket (20" X 20")	3522, 3523, 3525, 3533, 3622, 3623, 3633
 <b>CSG032</b>	536651	Door Gasket (24" X 36")	MC3600
CSK018	517295	Solenoid Valve Repair Kit	3523, 3533
CSK019	61301601738 & 518018 (old)	Dual Solenoid Valve Repair Kit	3522, 3533, 3622, 3623, 3633
CSK020	517298	Solenoid Valve Repair Kit	3523
CSK022	517297	Solenoid Valve Repair Kit	3523
CSK023	61301600442	Solenoid Valve Repair Kit	3522, 3523, 3525, 3533, 4200MPS, 6900AMPS
 <b>CSK024</b>	61301604084	Solenoid Valve Repair Kit	3622, 3623, 3633
CSK025	54186	Solenoid Valve Repair Kit	3622, 3623, 3633, 4200MPS, 6900AMPS
CSK048	68164	Solenoid Valve Repair Kit	3522, 3525, 3533
CSK050	54187	Solenoid Valve Repair Kit	3525, 4200MPS, 6900AMPS
CSK060	27557	Steam Trap Rebuild Kit	200, 3622, 3633, 4200MPS, 6900AMPS
CSK061	(No OEM Part # Available)	Solenoid Valve Repair Kit	3522, 3533, 3622, 3623, 3633
CSK062	(No OEM Part # Available)	Solenoid Valve Repair Kit	3523
CSK065	513700	Check Valve Rebuild Kit	3522, 3523, 3525, 3533, 3622, 3633
CSK066	513369	Check Valve Rebuild Kit	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633, 4200MPS, 6900AMPS
CSK067	(No OEM Part # Available)	Steam Trap Rebuild Kit	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633

RPI PART #	OEM PART #	DESCRIPTION	MODELS
CSP036	61301602085	PM Pkg (Solenoid Valves)	3522, 3523, 3533
CSP037	61301602090	PM Pkg (Solenoid Valves)	3525
CSP038	61301602091	PM Pkg (Filter & Screen)	3522, 3523, 3525, 3533
CSP039	61301602087	PM Pkg (Stm Trap & Chk Vlv)	3525
CSP040	61301602092	PM Pkg (Stm Trap & Chk Vlv)	3525
CSP041	61301602510	Spare Parts Package	3525
CSP042	61301602498	Spare Parts Package	3522
CSP043	61301602499	Spare Parts Package	3523
CSP044	61301602500	Spare Parts Package	3533
CSP045	61301602501	Spare Parts Package	3522
CSP046	61301602502	Spare Parts Package	3523
CSP047	61301602503	Spare Parts Package	3533
CSP068	61301602079	PM Pkg (Stm Trap & Chk Vlv)	3622, 3633
CSP069	(No OEM Part # Available)	PM Pkg (Stm Trap & Chk Vlv)	3522, 3523, 3525, 3533
CSP070	(No OEM Part # Available)	PM Pkg (Stm Trap & Chk Vlv)	3622, 3633
CSR028	86407	Solid State Relay	3522, 3523, 3525, 3533, 3622, 3623, 3633
CSR029	86408	Solid State Relay	3522, 3523, 3525, 3533, 3622, 3623, 3633
CSS053	43907	Screen Filter	3522, 3523, 3525, 3533, 3622, 3633
CSS054	45690	Screen, Detergent	3522, 3523, 3525, 3533, 3622, 3633
CST059	61301602492	Steam Trap	100, 200, 400/500, 700, 800, 3522, 3523, 3525, 3533, 3622, 3633
CSV004	517339	Solenoid Valve Assembly	3522, 3523, 3525, 3533, 3622, 3623, 3633
CSV005	517345	Solenoid Valve Assembly	3523
CSV007	517347	Solenoid Valve Assembly	3523, 3533
CSV008	517356	Dual Coil Solenoid Valve	3522, 3533, 3622, 3623, 3633
CSV009	517354	Solenoid Valve Assembly	3523
CSV010	517349	Solenoid Valve Assembly	3522, 3525
CSV011	517353	Solenoid Valve Assembly	3523
CSV012	541582	Solenoid Valve Assembly	3522, 3523, 3525, 3533, 4200MPS, 6900AMPS
CSV014	517342	Solenoid Valve Assembly	3622, 3623, 3633, 4200MPS, 6900AMPS
CSV033	56913	Safety Relief Valve	3622, 3633
CSV034	61301604180	Safety Relief Valve	3522, 3523, 3525, 3533
CSV035	532554	Relief Valve	3622, 3633
CSV056	14010	Check Valve	3522, 3523, 3525, 3533
CSV057	19114	Check Valve	3522, 3523, 3525, 3533
CSV058	08045	Check Valve	3522, 3523, 3525, 3533



# replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD  
U.S. POSTAGE PAID  
VAN NUYS, CA  
PERMIT NO. 2399

ADDRESS SERVICES REQUESTED

**LOOK INSIDE!**  
New Parts to fit  
CASTLE/GETINGE  
Bulk Sterilizers!