



CLEAR FILTER HOUSING NOW AVAILABLE TO FIT A-DEC UTILITY MODULE & REGULATORS!



RPI Part #ADK201
OEM Part #24-0228-00/
24-0234-00/030-019-00

FILTER KIT

Fits: Utility Module
(OEM Part #24-0352-00)
and Filter Regulators
(OEM Part #'s 24-0168-00,
24-0170-00 & 24-0175-00)
Models: All Models

The RPI Filter Housing is
CLEAR

so you can see when the filter
needs replacement without
having to remove the housing!

ADH200

ADF037

RP0325

All parts in Filter Kit also sold separately:


- Filter Housing (RPI Part #ADH200)
- Filter Element (RPI Part #ADF037)
- O-ring (RPI Part #RP0325)

UPGRADED RPI REPAIR KITS NOW INCLUDE BLACK AND WHITE VALVE BODIES!

RPI Part #ADK153
OEM Part # (No OEM Part # Available)

WHITE VALVE BODY REPAIR KIT

- This kit has been upgraded to include the White Valve Body (RPI Part #ADB197)
- Fits: Air Regulator Valve "White" Bodies
(OEM Part #24-0128-00 & #24-0135-00)
Models: Cascade, Decade, Excellence & Mini-Trol




- All parts in Repair Kit also sold separately:
- White Valve Body (RPI Part #ADB197)
 - Diaphragm (x2) (RPI Part #ADD019)
 - 9-Hole Gasket (RPI Part #ADG028)
 - 5-Hole Gasket (RPI Part #ADG027)
 - Conical Compression Spring (RPI Part #ADS102)
 - Plastic Piston (RPI Part #ADP157) with O-ring (RPI Part #RP0311)

RPI Part #ADK173
OEM Part # (No OEM Part # Available)

BLACK VALVE BODY REPAIR KIT

- This kit has been upgraded to include the Black Valve Body (RPI Part ADB198)
- Fits: Water Regulator Valve "Black" Bodies
(OEM Part #24-0129-00 & #24-0138-00)
Models: Cascade, Decade, Excellence & Mini-Trol



- All parts in Repair Kit also sold separately:
- Black Valve Body (RPI Part #ADB198)
 - Diaphragm (x2) (RPI Part #ADD019)
 - 9-Hole Gasket (RPI Part #ADG028)
 - Conical Compression Spring (RPI Part #ADS156)
 - Plastic Piston (RPI Part #ADP157) with O-ring (RPI Part #RP0311)