


# Control PC Boards to fit Old & New Style Midmark® • Ritter® UltraClave® M9 & M11 sterilizers!

Fits New Style M9 & M11



RPI Part #MIB129  
OEM Part #002-1992-00


**CONTROL PC BOARD (REFURBISHED)**

- Includes: Mounting Bracket
- Please Note: When shipping RPI your non-operational PC Board, **you must** include the Mounting Bracket in order to receive full credit!

**Fits:** Chassis

**Models:** M9 -020 & -022 Serial # V1697650 thru Present; and M11 -020 & -022 Serial # V1697650 thru Present

Fits New Style M9 & M11



RPI Part #MIB129  
OEM Part #022-0762-00


**CONTROL PC BOARD (REFURBISHED)**

- Includes: Mounting Bracket
- Please Note: When shipping RPI your non-operational PC Board, **you must** include the Mounting Bracket in order to receive full credit!

**Fits:** Chassis

**Models:** M9 -020 & -022 Serial # prefixes RN, RR & Serial # V2200 thru V1697649; and M11 -020 & -022 Serial # prefixes RS, RV & Serial # V2200 thru V1697649

Fits Old Style M9



RPI Part #MIB130  
OEM Part #002-0434-00/  
002-0434-03


**CONTROL PC BOARD (REFURBISHED)**

- Mounting Bracket not included
- Please Note: This PC Board **Does Not** include the Mounting Bracket. When shipping your non-operation PC Board to RPI, do not include your Mounting Bracket. You will need your Mounting Bracket to install the refurbished PC Board!

**Fits:** Chassis

**Models:** M9 -001 thru -019 with Serial # prefixes CZ, DA, DB, DX, DY, FD & OM

Fits Old Style M11



RPI Part #MIB131  
OEM Part #002-0501-00

**CONTROL PC BOARD (REFURBISHED)**

- Mounting Bracket not included
- Please Note: This PC Board **Does Not** include the Mounting Bracket. When shipping your non-operation PC Board to RPI, do not include your Mounting Bracket. You will need your Mounting Bracket to install the refurbished PC Board!

**Fits:** Chassis

**Models:** M11 -001 thru -019 with Serial # prefixes ES, ET, FP, FR & GB

## About RPI's Control PC Board Exchange Program

**Here's How the \$100 Credit Exchange Program Works.**

- 1) Place an order for a Control PC Board (PCB), and RPI will send a shipping box and prepaid UPS return label to you. You will not pay for shipping when using this box.
- 2) Using the shipping box and UPS return label, ship as follows:
  - **Include** the PCB and Mounting Bracket if shipping the *New Style* M9/M11 PCB (OEM Part #002-1992-00/022-0762-00). If you do not include the Mounting Bracket, you will only receive a \$75 Credit instead of \$100.
  - **Do not include** the Mounting Bracket if shipping the *Old Style* M9/M11 PCB (OEM Part #'s 002-0434-00/002-0434-03/002-0501-00).
- 3) Once the non-operational board has been received by RPI and verified that it is in reasonable condition without burn holes, a replacement PCB will be shipped to you. Your invoice will reflect the cost of the board plus the \$100 credit for the Core.

**Important Information.** PC Board performance verified and tested and Quality Control test printout is included with each PC Board. 90 Day Warranty.

**Warning.** To avoid damaging the replacement PC Board and possibly voiding the warranty, check the condition of the *Heater and/or Pulse Solenoid* prior to installing the replacement PC Board. DO NOT use PC Board in 220V or version "D" sterilizers.




**You asked for It ... RPI Now Offers the Red Colored PC Board to fit newer Midmark • Ritter UltraClave M9 & M11 sterilizers with a Green LCD Display!**

RPI now offers the **red colored PC Boards** that fit the newer Midmark sterilizers with a green display through our PC Board Exchange Program - receive a \$100\* credit with the bracket, and \$75.00\* without the bracket. These boards fit the M9 -020 thru -022 Serial # V1697650 thru Present and the M11 -020 thru -022 Serial # V1697650 thru Present. Just ask for the Control PC Board (RPI Part #MIB132).

Don't forget, we already have a program for the **green colored PC Boards** that also fit newer Midmark sterilizers with a green display. This PCB fits the M9 -020 thru -022 with serial number prefixes RN, RR, and Serial # V2200 thru V1697649; and the M11 -020 thru -022 with serial number prefixes RS, RV, and Serial # V2200 thru V1697649. Just ask for RPI Part #MIB129. Also note that both the red and green boards will only be exchanged for the exact same style board - red for red, and green for green - this is the only way in which we can keep this program working.

\* Credit is given after RPI has inspected the PC Board and found it to be in good condition.



**RPI Part #MIB132**

**PC Boards to fit Midmark® • Ritter® UltraClave® M9 & M11 Sterilizers!**


ADDRESS SERVICES REQUESTED

**replacement parts industries, inc.**  
 625 Cochran Street • Simi Valley, California • 93065




PERMITTED STANDARD  
 U.S. POSTAGE PAID  
 VAN NUYS, CA  
 PERMIT NO. 2399


**RPI now offers four Control PC Boards to fit Midmark® • Ritter® UltraClave® M9 & M11 Sterilizers!**




**RPI Part #MIB132**  
 Fits new style  
**M9 & M11**



**RPI Part #MIB129**  
 Fits new style  
**M9 & M11**



**RPI Part #MIB130**  
 Fits old style  
**M9**



**RPI Part #MIB131**  
 Fits old style  
**M11**