



**Parts to fit
Sterrad NX, 100NX & 100S
Hydrogen Peroxide Sterilization Systems!**



**RPI Part #SDK100
OEM Part #35-53016-001**

NEUTRAL DENSITY FILTER KIT

- Includes all Neutral Density Filters necessary to complete the Optical Density Calibration procedure
- Includes individual certification of calibration for each lens
- Includes all parts as shown
- Also available: Optical Window Kit (RPI Part #SDK014), Detector Window Kit (RPI Part #SDK116), and UV Lamp Assembly (RPI Part #SDA006)

Fits: Optical Detector Assembly
Models: NX & 100NX



Neutral Density Filter Kit (RPI Part #SDK100) Includes All Parts Shown Below!



**RPI Part #SDF101
OEM Part # (No OEM Part # Available)**

NEUTRAL DENSITY FILTER (#1)

- 0.0400 Optical Density
- Fits:** Optical Detector Assembly
Models: NX & 100NX



**RPI Part #SDF103
OEM Part # (No OEM Part # Available)**

NEUTRAL DENSITY FILTER (#3)

- 0.2285 Optical Density
- Fits:** Optical Detector Assembly
Models: NX & 100NX



**RPI Part #SDF102
OEM Part # (No OEM Part # Available)**

NEUTRAL DENSITY FILTER (#2)

- 0.1475 Optical Density
- Fits:** Optical Detector Assembly
Models: NX & 100NX



**RPI Part #SDF104
OEM Part # (No OEM Part # Available)**

NEUTRAL DENSITY FILTER (#4)

- 0.4020 Optical Density
- Fits:** Optical Detector Assembly
Models: NX & 100NX



**RPI Part #RXT109
OEM Part # (No OEM Part # Available)**

OPTICAL LENS TOWLETTE

- 4 per package
 - Anti-static, anti-fog and Anti-streaking single use premoistened wipe
- Models:** NX & 100NX



**RPI Part #RXC110
OEM Part # (No OEM Part # Available)**

STORAGE CASE

- Size: 7.91" x 4.25" x 2.04"
- Models:** NX & 100NX



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 2002SD Rd5 SHOW 07/21

New Parts to fit Sterrad NX & 100NX H2O2 Sterilization Systems



RPI Part #SDK116
OEM Part #32-03557-0-001

DETECTOR WINDOW KIT

- Includes: Lens and Viton® (FKM) O-Ring
- Fits:** Optical Detector Assembly
Models: NX & 100NX

LOOK INSIDE!

Also available: *Neutral Density Filter Kit*
 (RPI Part #SDK100) and *More New Parts!*



RPI Part #SDK014
OEM Part #32-51413-0-001/22-00537-209

OPTICAL WINDOW KIT

- Includes: (2) Optical windows and (1) Viton® (FKM) O-Ring
 - Included in: PM-2 Kit (RPI Part #SDK003)
- Fits:** Hydrogen Peroxide Monitor and Optical Detector Assembly
Models: NX & 100NX



RPI Part #SDA006
OEM Part #40-52125-001

UV LAMP ASSEMBLY

- Rated at 900 VRMS
 - 12" lg. cable with connector
 - Factory calibrated
 - Included in: PM-2 Kit (RPI Part #SDK003)
- Fits:** Hydrogen Peroxide Monitor and Optical Detector Assembly
Models: NX & 100NX

New Parts to fit Sterrad NX & 100NX H2O2 Sterilization Systems


ADDRESS SERVICES REQUESTED

PERMITTED NO. 2399
 VAN NUYS, CA
 U.S. POSTAGE PAID
 PRESORTED STANDARD

625 Cochran Street • Simi Valley, California • 93065
replacement parts industries, inc.



Also Available to fit the 100S




RPI Part #SDK117
OEM Part # (No OEM Part # Available)

SERVICE ACCESS KEY

- Used with Auto Test to gain access to the Service Mode
- Fits:** J2 Harness Connector
OEM: STERRAD
Models: 100S

Don't Forget ...

Essential tool to protect sensitive equipment during maintenance!



RPI Part #RPS998
OEM Part # (No OEM Part # Available)

GROUNDING STRAP

- For grounding to an ESD work bench or a chassis
 - Includes: Strap, 6' long coiled cord with banana plug and detachable alligator clamp
- Models:** NX & 100NX