




New line of parts to fit Welch Allyn Wall Transformers!



RPI Part #WAB001
OEM Part #747284-501
POWER SUPPLY BOARD (3.5V)

- 3.5VDC
- Includes all parts as shown
- Also available: Position Control Knob (RPI Part #BUK043)



Fits: Chassis
Model: 74710



RPI Part #WAK003
OEM Part #747300-501 (Complete Handpiece Assembly)
COILED CORD KIT

- Size: 2 Conductor #20 AWG x ~3 ft. lg.
- Includes all parts as shown
- Now you can replace just the Coiled Cord instead of the whole Handpiece Assembly - It's the RPI Advantage!
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)


Fits: Handpiece
Model: 74710



RPI Part #WAB004
OEM Part #747145
WALL BRACKET

- Size: 4.16" lg. x 3.0" wd. x .06" thk.
- Material: Aluminum
- Includes: Mounting hardware

Fits: Mounting Plate
Model: 74710



RPI Part #WAB007
OEM Part #767035-501
WALL BRACKET

- Size: 6.0" lg. x 2.2" wd. x .06" thk.
- Material: Aluminum
- Includes: Mounting hardware

Fits: Mounting Plate
Model: 767




RPI Part #BUK043
OEM Part #747116-1



POSITION CONTROL KNOB

- Material: Black plastic knob, aluminum face plate
- Size: .92" OD x .64" lg.
- Shaft Size: 1/4" dia.
- Includes: (2) installed set screws

Fits: Control Panel
Model: 74710




RPI Part #WAS002
OEM Part #747014-2



MICROSWITCH

- Momentary
- 4A @ 250VAC; SPDT
- .058" dia. Pin Terminals

Fits: Handle Mounting Bracket
Model: 74710




RPI Part #TUS009
OEM Part #483311-2



POWER SWITCH

- Rocker switch (On-Off)
- 20A @ 125/250VAC; DPST
- Color: Black body with illuminated green lens; Legend: I/O
- .25" Tabs

Fits: Main Housing
Model: 74710




RPI Part #WAS008
OEM Part #236707-1



POWER SWITCH

- Rocker Switch (ON - OFF)
- 16A @ 250VAC; DPST
- Color: Black body with illuminated green lens
- .25" Tabs

Fits: Main Housing
Model: 767



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Welch Allyn equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. Show 09/14




replacement parts industries, inc.

20338 Corisco Street • Chatsworth, California • 91311

PRESORTED STANDARD
U.S. POSTAGE PAID
VAN NUYS, CA
PERMIT NO. 2399

ADDRESS SERVICES REQUESTED

RPI Offers New Line of Parts to fit Welch Allyn Wall Transformers!



RPI Part #CAK120
OEM Part #747306-501 (Power Cord only)/711314 (Strain Relief only)

POWER CORD KIT

- Length: 8 ft.
- Color: Gray
- 18/3 SJT, 60°C; 10A @ 125VAC
- Includes all parts as shown
- Also available: Strain Relief Bushing Pliers (RPI Part #RPT172)


Fits: Chassis
Model: 74710

RPI Part #RPC557
OEM Part #761069

HOSPITAL GRADE POWER CORD

- Length: 8 ft.; Color: Gray
- 16/3 SJT, 60°C; 13A @ 125VAC
- Plug: NEMA 5-15P
- Connector: EN60320-C13

Fits: Power Input Module
Model: 767




RPI Part #RPF508
OEM Part #236706-3111

FUSE (1/4A, 250V) - FAST ACTING

- 5 per package
- Type: Glass
- Size: 5mm dia. x 20mm lg.
- Included in: Fuse Kit (RPI Part #RPK520)

Fits: Power Supply Board
Model: 767



Lamps for Exam Heads


RPI Part #LMP109
OEM Part #WA-03100-U

LAMP - 6/pkg

- 3.5V, 2.52W
- Shape: T1.5; Xenon

Fits: Examination Head

Models: 20000 Otoscope, 20200 Otoscope, 20202 Otoscope, 21700 Otoscope, 25000 Otoscope, 25200 Otoscope, 25020 Otoscope, 73500 Handle Adapter, 28100 Tongue Blade Holder, 26530 Illuminator, 27000 Illuminator, 27050 Illuminator, 41100 Illuminator & 43300 Illuminator




RPI Part #LMP110
OEM Part #WA-03000-U

LAMP - 6/pkg

- 3.5V, 2.52W
- Shape: T1.5; Xenon
- Single contact 8mm flange

Fits: Examination Head

Models: 11600 Ophthalmoscope, 11605 Ophthalmoscope, 11610 Ophthalmoscope, 11710 Ophthalmoscope, 12400 Strabismoscope, 43700 Episcopes & 18000 Retinoscope




RPI Part #LMP111
OEM Part #WA-06500-U

LAMP - 6/pkg

- 3.5V, 2.73W
- Shape: TL1.5; Xenon

Fits: Examination Head

Models: 23810 Otoscope & 23820 Otoscope



RPI Part #LMP112
OEM Part #WA-03800-U

LAMP - 6/pkg

- 3.5V, 3.325W
- Shape: T1.5; Xenon

Fits: Examination Head

Models: 11800 Ophthalmoscope, 11810 Ophthalmoscope & 11820 Ophthalmoscope

