



CABLES TO FIT GE ACQUISITION MODULES

RPI PART #	DESCRIPTION	LENGTH
<p>GEC020</p>	<p>ACQUISITION MODULE CABLE</p> <ul style="list-style-type: none"> • OEM Part #700044-201, 700044-202 & 700044-204 • Fits: AM4 Acquisition Module & AM5 Acquisition Module 	<p>~17 ft. overall length</p>  <p>The diagram shows a cable with a multi-pin connector on the left, a central section of twisted wires, a long section of braided shielded cable, and a multi-pin connector on the right.</p>
<p>GEC021</p>	<p>ACQUISITION MODULE CABLE</p> <ul style="list-style-type: none"> • OEM Part #2016560-002, 2016560-003, 700657-002 & 700657-003 • Fits: CAM 14 Acquisition Module 	<p>~17 ft. overall length</p>  <p>The diagram shows a cable with a circular multi-pin connector on the left, a central section of twisted wires, a long section of braided shielded cable, and a multi-pin connector on the right.</p>