
























DIN/UNSHIELDED PATIENT CABLES

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	SUGGESTED LEADWIRE KITS
KCA001	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	6-Pin; 1K Ω Resistors Jumper E to C	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA004	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	6-Pin; 1K Ω Resistors	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA007	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	6-PIN Jumper E to B	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA010	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	12-Pin; 1K Ω Resistors	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA013	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	5-Pin	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM008 & LKM011 Disposable Leadwire Sets RPI Part #'s LDM302 & LDM309 Veterinary Leadwire Set RPI Part #LVM702
KCA016	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	6-Pin	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA019	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	11-Pin jumper 5 to 6 with 120M Ω Resistor	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA022	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	8-Pin; 1K Ω Resistors	Din/Unshielded .060" Pin (Safety) 	Reusable Leadwire Sets RPI Part #'s LKM008 & LKM011 Disposable Leadwire Sets RPI Part #'s LDM302 & LDM309 Veterinary Leadwire Set RPI Part #LVM702 NOTES: <ul style="list-style-type: none"> This cable replaces RPI Part #HPC004

DIN/UNSHIELDED PATIENT CABLES *(Continued)*

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END	SUGGESTED LEADWIRE KITS
KCA025	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	8-Pin; 1K Ω Resistors 	Din/Unshielded .060" Pin (Safety)  WHITE RED BLACK 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701
KCA035	PATIENT CABLE - 3 LEAD DIN	3	~10 ft.	12-Pin Right Angle; 10K Ω Resistors 	Din/Unshielded .060" Pin (Safety)  WHITE RED BLACK 	Reusable Leadwire Sets RPI Part #'s LKM002 & LKM005 Disposable Leadwire Sets RPI Part #'s LDM301 & LDM308 Veterinary Leadwire Set RPI Part #LVM701 NOTES: • OEM Part #11110-000029 & 3006218-02
KEA001	PATIENT CABLE - 5 LEAD DIN	5	~10 ft.	11-Pin 	Din/Unshielded .060" Pin (Safety)  WHITE GREEN BROWN RED BLACK 	Reusable Leadwire Sets RPI Part #'s LKM026, LKM029 & LKM030 Disposable Leadwire Sets RPI Part #'s LDM303 & LDM310 Veterinary Leadwire Set RPI Part #LVM703
KEA005	PATIENT CABLE - 5 LEAD DIN	5	~10 ft.	12-Pin; 1K Ω Resistors 	Din/Unshielded .060" Pin (Safety)  WHITE GREEN BROWN RED BLACK 	Reusable Leadwire Sets RPI Part #'s LKM026, LKM029 & LKM030 Disposable Leadwire Sets RPI Part #'s LDM303 & LDM310 Veterinary Leadwire Set RPI Part #LVM703
KEA017	PATIENT CABLE - 5 LEAD DIN	5	~10 ft.	8-Pin; 1K Ω Resistors 	Din/Unshielded .060" Pin (Safety)  WHITE GREEN BROWN RED BLACK 	Reusable Leadwire Sets RPI Part #'s LKM026, LKM029 & LKM030 Disposable Leadwire Sets RPI Part #'s LDM303 & LDM310 Veterinary Leadwire Set RPI Part #LVM703 NOTES: • This cable replaces RPI Part #HPC006