

LISTING OF RPI PARTS THAT FIT TUTTNAUER STERILIZERS			1730E	1730EK	1730M	1730MK	2340E	2340EK	2340M	2340MK	2540E	2540EK	2540M	2540MK	3870E	3870EK	3870M	3870MK
OEM PART #	RPI PART #	DESCRIPTION																
01620004	TUT038	SAFETY THERMOSTAT (MANUAL RESET)			•	•			•	•			•	•				
01620101	TUT007	CONTROL THERMOSTAT			(8)				(8)				(8)					
01620103	TUT039	SAFETY THERMOSTAT (AUTOMATIC RESET)			•	•			•	•			•	•				
01620301	TUS035	OVER TEMP SAFETY SWITCH		(7)		(7)		(7)		(7)		(7)		(7)				
01720001	TUH027	HEATER ELEMENT	•		•													
01720002	TUH005	HEATER ELEMENT					•		•									
01720003	TUH006	HEATER ELEMENT									•		•					
01720011	TUH004	HEATER ELEMENT		(2)		(2)												
01720012	TUH016	HEATER ELEMENT		(3)		(3)												
01720013	TUH017	HEATER ELEMENT						(3)		(3)								
01720014	TUH018	HEATER ELEMENT										(3)		(3)				
01910011	TUT008	TIMER			(6)	(6)			(6)	(6)			(6)	(6)				
01910097	TUB024	CIRCUIT BREAKER	(1)	(3)	(1)	(3)	(3)		(3)		(3)		(3)					
01910098	TUB023	CIRCUIT BREAKER		(2)		(2)	(2)	(3)	(2)	(3)	(2)	(3)	(2)	(3)	(3)			(3)
01910172	TUS009	POWER SWITCH			•	•			•	•			•	•				
01910190	TUS014	DOOR SWITCH	•	•	•	•	•	•	•	•	•	•	•	•				
01910191	TUS013	MICROSWITCH			•	•			•	•			•	•				
01910257	MDL021	SIGNAL LIGHT - ORANGE			(2)				(2)				(2)					
01910258	MDL022	SIGNAL LIGHT - GREEN			(2)				(2)				(2)					
01960007	MDL019	SIGNAL LIGHT - RED			(2)				(2)				(2)					
02300011	TUG012	PRESSURE GAUGE WITH MAX INDICATING POINTER							•				•	•				
02610005	TUG002	DOOR GASKET							(4)									
02610020	TUG001	DOOR GASKET	•	•	•	•												
02610023	TUG003	DOOR GASKET									•	•	•					
02610118	TUG021	DOOR GASKET					•	•	(5)	(5)								
02819995	TUC028	POWER CORD	•	•	•	•	•	•	•	•	•	•	•	•				
03110002	TUV011	SAFETY RELEASE VALVE	•	•	•	•	•	•	•	•	•	•	•	•				
CB842010	TUJ034	AIR VALVE JET		•		•		•		•		•		•		•		•
CT241111	TUK030	DOOR BELLOWS ASSEMBLY KIT	•	•	•	•	•	•	•	•	•	•	•	•				
CT842010	TUJ033	AIR VALVE JET	•		•		•		•		•		•		•		•	
N/A	TUB019	DOOR BELLOWS	•	•	•	•	•	•	•	•	•	•	•	•				
N/A	TUH015	ATTACHING HARDWARE	•	•	•	•	•	•	•	•	•	•	•	•				

(1) Fits 110V and 220V units

(2) Fits 110V units

(3) Fits 220V units

(4) Fits units with S/N 8805 and below

(5) Fits units with S/N 8806 and above

(6) Old OEM Part #01910005

(7) Only fits Models prior to 1993

(8) Only fits manual Models