



New Parts to fit CASTLE/GETINGE

Other Solenoid Valve Assemblies & Kits Available

SOL-1 STEAM TO CHAMBER

RPI Part #CSV004 OEM Part #517339
 Models: 3522, 3523, 3525, 3533, 3622, 3623 & 3633
 See also: Valve Repair Kit (RPI Part #AMK177)

SOL-1 STEAM TO CHAMBER

RPI Part #CSV005 OEM Part #517345
 Model: 3523
 See also: Valve Repair Kit (RPI Part #AMK183)

SOL-3 AIR-IN

RPI Part #CSV007 OEM Part #517347
 Models: 3523 & 3533
 See also: Valve Repair Kit (RPI Part #CSK018)

SOL-4 MULTI BYPASS

RPI Part #CSV004 OEM Part #517339
 Models: 3523, 3622, 3623 & 3633
 See also: Valve Repair Kit (RPI Part #AMK177)

SOL-5 SLOW EXHAUST

RPI Part #CSV009 OEM Part #517354
 Model: 3523
 See also: Valve Repair Kit (RPI Part #CSK020)

SOL-5 & 6 SLOW EXHAUST & VAC-WATER DUAL SOLENOID VALVE REPAIR KIT

RPI Part #CSK019 OEM Part #61301601738
 Models: 3522, 3533, 3622, 3623 & 3633

SOL-6 VAC-WATER

RPI Part #CSV010 OEM Part #517349
 Models: 3522 & 3525
 See also: Valve Repair Kit (RPI Part #AMK174)

SOL-7 MULTI BYPASS SOLENOID VALVE REPAIR KIT

RPI Part #CSK022 OEM Part #517297
 Model: 3523

SOL-7 MULTI STEAM

RPI Part #CSV004 OEM Part #517339
 Models: 3622, 3623 & 3633
 See also: Valve Repair Kit (RPI Part #AMK177)

SOL-8 JACKET


RPI Part #CSV004 OEM Part #517339
 Models: 3522, 3523, 3533, 3622, 3623 & 3633
 See also: Valve Repair Kit (RPI Part #AMK177)

SOL-8 JACKET

RPI Part #CSV014 OEM Part #517342
 Models: 3622, 3623, 3633, 4200 & 6900
 See also: Valve Repair Kit (RPI Part #CSK025)

SOL-10 & 11 DOOR

RPI Part #CSV012 OEM Part #541582
 Models: 3522, 3523, 3525, 3533, 4200 & 6900
 See also: Valve Repair Kit (RPI Part #CSK023)




SOL-5 & 6 Slow Exhaust & Vac-Water

RPI Part #CSV008
 OEM Part #517356

DUAL COIL SOLENOID VALVE

- This is **not** an ASCO® Valve Assembly
 - Includes all parts as shown
- Models: 3522, 3533, 3622, 3623 & 3633




RPI Part #CSK061
 OEM Part #

(No OEM Part # Available)

SOLENOID VALVE REPAIR KIT

- Fits RPI Dual Coil Solenoid Valve (RPI Part #CSV008) **only**; For kit to repair OEM Dual Coil Solenoid Valve NO. 517356, see RPI Part #CSK019
 - Includes all parts as shown
- Models: 3522, 3533, 3622, 3623 & 3633


Also available, **DUAL SOLENOID VALVE REPAIR KIT** (RPI Part #CSK019) to fit Castle/Getinge valve #517356 **only**.



RPI Part #CSC063
 OEM Part # (No OEM Part # Available)

COIL

- Fits Dual Coil Solenoid Valve (RPI Part #CSV008) **only**
- Models: 3522, 3533, 3622, 3623 & 3633




SOL-7 Multitherm Steam

RPI Part #CSV011
 OEM Part #517353

SOLENOID VALVE ASSEMBLY

- This is **not** an ASCO® Valve Assembly
 - Includes all parts as shown
- Model: 3523 (Sol-7 Multitherm Steam)




RPI Part #CSK062
 OEM Part #

(No OEM Part # Available)

SOLENOID VALVE REPAIR KIT

- Fits Solenoid Valve Assembly (RPI Part #CSV011) **only**
 - Includes all parts as shown
- Model: 3523 (Sol-7 Multitherm Steam)



RPI Part #CSC064
 OEM Part # (No OEM Part # Available)

COIL

- Fits Solenoid Valve Assembly (RPI Part #CSV011) **only**
- Model: 3523 (Sol-7 Multitherm Steam)



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Castle/Getinge equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB 03/09