



... Proud to be a **Healthcare Components Group** company ...



**HEALTHCARE
Components Group**

MIDMARK • RITTER



New Autoclave Data Loggers to fit Midmark® • Ritter® M9 & M11

Automatic • Wireless • Paperless

Securely store autoclave cycle data to the cloud and on a removable SD card and eliminate the need for printing.

Features

- Wireless Data Recorder with LCD Display
- Removeable 16gb micro SD card
- Secure Wi-Fi
- Password protected
- FCC Certified
- Unique naming of each autoclave
- 5 year time and date battery
- Automatic time keeping and daylight savings time adjustment
- Wireless upload to offsite Dropbox™ cloud storage



**RPI Part #MIL337
OEM Part #9A682001**



DATA LOGGER (M19)

- The Autoclave Logger Thing M19
- Size: 2-1/2" x 4" x 1"
- Includes: Removable 16gb micro SD card, SD card adaptor, stand, cable, and user/setup manual.

Fits: Autoclave

Models: M9-001 thru -019 & M11-001 thru -019



**RPI Part #MIL338
OEM Part #9A682001**



DATA LOGGER (M22)

- Autoclave Logger Thing M22
- Size: 2-1/2" x 4" x 1"
- Includes: Removable 16gb micro SD card, SD card adaptor, stand, cable, user/setup manual, and drill template for use with printer cover plate.

Fits: Autoclave

Models: M9-020 thru -022 & M11-020 thru -022



**RPI Part #MIL340
OEM Part #9A682001**



DATA LOGGER (M42)

- Autoclave Logger Thing M42
- Size: 2-1/2" x 4" x 1"
- Includes: Removable 16gb micro SD card, SD card adaptor, stand, cable, user/setup manual, and drill template for use with printer cover plate.

Fits: Autoclave

Models: M9-040 thru -043 & M11-040 thru -043



replacement parts industries, inc. A Healthcare Components Group Company

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. 2104DL Web 01/23

SECTION

10

PAGE

1026

UPDATED

01/23