



PM Packs & More

to fit Getinge/Castle® 400 & 500 Series Steam Sterilizers



RPI Part #CSP110
OEM Part #61301605861

PM PACK (LOWER)


- **PM Packs include rebuild kits for all valves.**
 - Includes: Solenoid Valve Repair Kit (RPI Part #CSK090), Solenoid Valve Repair Kit (RPI Part #CSK019), (4) Check Valve Rebuild Kit (RPI Part #CSK066), Solenoid Valve Repair Kit (RPI Part #CSK096), (2) Solenoid Valve Repair Kit (RPI Part #CSK094), Solenoid Valve Repair Kit (RPI Part #AMK177), Steam Trap (RPI Part #CST059), and Strainer Element (RPI Part #CSS114)
 - Also available: PM Pack (Upper) (RPI Part #CSP111)
- Models:** 500 Series



RPI Part #CSP112
OEM Part #61301605863

PM PACK (LOWER)

- **PM Packs include rebuild kits for all valves.**
 - Includes: Solenoid Valve Repair Kit (RPI Part #CSK090), Solenoid Valve Repair Kit (RPI Part #CSK019), (5) Check Valve Rebuild Kit (RPI Part #CSK066), Solenoid Valve Repair Kit (RPI Part #CSK092), (2) Solenoid Valve Repair Kit (RPI Part #CSK094), Solenoid Valve Repair Kit (RPI Part #AMK177), Steam Trap (RPI Part #CST059), and Strainer Element (RPI Part #CSS114)
 - Also available: PM Pack (Upper) (RPI Part #CSP111)
- Models:** 400 Series



RPI Part #CSP111
OEM Part #61301605862

PM PACK (UPPER)

- **PM Packs include rebuild kits for all valves.**
 - Includes: (3) Solenoid Valve Repair Kit (RPI Part #AMK177), Solenoid Valve Repair Kit (RPI Part #CSK108), (2) Check Valve Rebuild Kit (RPI Part #CSK066), Filter, Air-In (RPI Part #CSF052), Steam Trap (RPI Part #CST082), Strainer Element (RPI Part #CSS114), and (2) Safety Valve (45 PSI) (RPI Part #CSV084)
 - Also available: PM Pack (Lower) (RPI Part #CSP110) and PM Pack (Lower) (RPI Part #CSP112)
- Models:** 400 Series & 500 Series



RPI Part #CSP113
OEM Part #61301605864

PM PACK (DOUBLE DOOR)

- **PM Packs include rebuild kits for all valves.**
 - Includes: (2) Solenoid Valve Repair Kit (RPI Part #AMK177) and (2) Check Valve Rebuild Kit (RPI Part #CSK066)
 - Also available: PM Pack (Lower) (RPI Part #CSP110), PM Pack (Upper) (RPI Part #CSP111), and PM Pack (Lower) (RPI Part #CSP112)
- Models:** 400 Series & 500 Series





RPI Part #CSV089
OEM Part #61301601701

SOLENOID VALVE

- Ports: 3/4" FPT, 2-Way: N.C.
 - 24VDC 60 Hz./22VDC 50 Hz., 10.1 W (3.6 Ohms)
 - Seal Material: PTFE and EPDM
 - Also available: Solenoid Valve Repair Kit (RPI Part #CSK090)
- Fits:** Sol-3 Fast Exhaust & Sol-2 Chamber Drain
Models: 400 Series & 500 Series

Repair Kit



RPI Part #CSK090
OEM Part #61301600443/764337708G

SOLENOID VALVE REPAIR KIT

- Seal Material: PTFE and EPDM
 - Includes all parts as shown
 - Also available: Solenoid Valve (RPI Part #CSV089)
- Fits:** Sol-3 Fast Exhaust & Sol-2 Chamber Drain
Models: 400 Series & 500 Series



RPI Part #CSV091
OEM Part #61301601696

SOLENOID VALVE

- Ports: 3/4" FPT, N.C.
 - 24VDC 60 Hz./22VDC 50 Hz., 6.1 W (8.3 Ohms)
 - Seal Material: PTFE and EPDM
 - Also available: Dual Solenoid Valve Repair Kit (RPI Part #CSK019)
- Fits:** Sol-5/Sol-6 Ejector and Cooling
Models: 400 Series & 500 Series

Repair Kit



RPI Part #CSK019
OEM Part #518018/61301601738

DUAL SOLENOID VALVE REPAIR KIT

- Repair kit for use with Sol-5 & Sol-6 (Getinge/Castle Valve No. 517356) only
 - For **Green** Coil Valves
 - Includes: Solenoid "A": Diaphragm style repair kit, EPDM seal and Solenoid "B": Piston style valve, PTFE seal
 - Also available: Solenoid Valve (RPI Part #CSV091)
- Fits:** Sol-5/Sol-6 Ejector and Cooling
Models: 100 Series, 200 Series, 400 Series, 500 Series, 700 Series, 3522, 3533, 3622, 3623 & 3633



RPI Part #CSS114
OEM Part #513448

STRAINER ELEMENT

- Material: Stainless steel
 - Size: 1.03" OD x 2.45" lg.
- Fits:** 3/4" Strainer
Models: 400 Series & 500 Series

RPI Part #CSF052
OEM Part #61301600303

FILTER, AIR-IN

- 2.84" OD x 3.10" lg.
 - 3/4" FPT
 - 0.03 micron HEPA media
- Fits:** Air Supply Line
Models: 100 Series, 200 Series, 3522, 3523, 3525, 3533, 3622, 3623, 3633, 400 Series, 4200MPS, 500 Series, 6900 Series, 700 Series, 800-42" & 800-76"





RPI Part #AMK177
OEM Part #61301601739/764265001

SOLENOID VALVE REPAIR KIT

- For 3/8" FPT Port Valves
 - PTFE seal
 - Includes all parts as shown
 - Also available: Solenoid Valve Assembly (RPI Part #CSV004)
- Fits:** Sol-1 Steam to Chamber, Sol-4 Multitherm Valve, Sol-7 Multitherm, Sol-8 Steam to Jacket, Sol-10 Door Seal/Unseal, Sol-12 Door Seal/Unseal, Sol-16 Vent
- Models:** 100 Series, 3522, 3523, 3525, 3533, 3622, 3623, 3633, 400 Series, 500 Series & 700 Series



RPI Part #CSK092
OEM Part # (No OEM Part # Available)

SOLENOID VALVE REPAIR KIT

- Fits valve's upper section and manual override body
 - Viton® seals
 - For **Black** Coil Valves
 - Includes all parts as shown
 - Now you can repair the valve with this one-of-a-kind RPI valve repair kit. It's the RPI Advantage!
 - Also available: Solenoid Valve Repair Kit (RPI Part #CSK096) - Fits bottom section
- Fits:** Sol-6 Water to Ejector (Black Coil) Valve Assembly (OEM Part #61301609315)
- Models:** 400 Series, 500 Series, 700 Series & 800-42"



RPI Part #CSK066
OEM Part #513369

CHECK VALVE REBUILD KIT

- For 1/2", 3/8" or 3/4" style valve
 - Material: Viton® O-ring; Stainless steel springs and retainer; Brass disc holder
 - Includes all parts as shown
 - This Kit includes two different seat and spring sets for two different applications
- Fits:** Check Valve (OEM Part #'s 513361, 513362 and 513363)
- Models:** 100 Series, 200 Series, 3522, 3523, 3525, 3533, 3622, 3633, 400 Series, 4200MPS, 500 Series, 6900 Series, 700 Series, 800-42" & 800-76"



RPI Part #CSK096
OEM Part #61301609505

SOLENOID VALVE REPAIR KIT

- Fits valve's lower section
 - Viton® seals
 - For **Black** Coil Valves
 - Includes: O-ring (RPI Part #RPO864), spring, diaphragm assembly and (2) O-rings
 - Also available: Solenoid Valve Repair Kit (RPI Part #CSK092) - Fits top section
- Fits:** Sol-6 Water to Ejector (Black Coil) Valve Assembly (OEM Part #61301609315)
- Models:** 400 Series, 500 Series, 700 Series & 800-42"



RPI Part #CSK094
OEM Part #61301608963

SOLENOID VALVE REPAIR KIT

- Teflon® seal
 - For **Black** Coil Valves
 - Includes all parts as shown
- Fits:** Sol-5 Cooling Water (Black Coil) Valve and Sol-14 Secondary Drain (Black Coil) Valve
- Models:** 400 Series, 500 Series, 700 Series & 800-42"



RPI Part #CSK108
OEM Part #61301609451

SOLENOID VALVE REPAIR KIT

- Viton® seals
 - For **Black** Coil Valves
 - Includes all parts as shown
- Fits:** Sol-3 Air-In (Black Coil) Valve Assembly
- Models:** 400 Series & 500 Series



RPI Part #CST059
OEM Part #61301602492

STEAM TRAP

- 1/2" MPT x 1/2" FPT
- Material: Brass body; stainless steel seat
- Max. rated pressure: 125 PSI;
Max. operating temperature: 450°F
- Also available: Steam Trap Rebuild Kit (RPI Part #CSK067)

Fits: Jacket or Chamber to Jacket

Models: 100 Series, 200 Series, 3522, 3523, 3525, 3533, 3622, 3633, 400 Series, 500 Series, 700 Series, 800-42" & 800-76"



RPI Part #CST082
OEM Part #68722

STEAM TRAP

- 3/8" FPT
- Material: Stainless steel body and disc
- Max operating pressure: 600 PSI;
Minimum operating pressure: 3.5 PSI
- Max allowable temperature: 800°F
- Includes all parts as shown

Fits: Steam Line

Models: 400 Series, 4200MPS, 500 Series & 700 Series



RPI Part #CSV084
OEM Part #61301608100

SOLENOID VALVE (45 PSI)

- 1/2" MPT x 1/2" FPT
- Material: Bronze body; EPDM seal
- Cracking pressure: 45 PSI

Fits: Chamber Jacket

Models: 400 Series, 500 Series & 700 Series

PM & Spare Parts packages for Getinge/Castle Bulk Sterilizers

RPI Part #	OEM Part #	Description	Fits Model(s)
CSP036	61301602085	PM PKG (SOLENOID VALVES)	3522, 3523 & 3533
CSP037	61301602090	PM PKG (SOLENOID VALVES)	3525
CSP038	61301602091	PM PKG (FILTER & SCREEN)	3522, 3523, 3525 & 3533
CSP039	61301602087	PM PKG (STM TRAP & CHK VLV)	3525
CSP040	61301602092	PM PKG (STM TRAP & CHK VLV)	3525
CSP068	61301602079	PM PKG (STM TRAP & CHK VLV)	3622 & 3633
CSP069	<i>No OEM Part # Available</i>	PM PKG (STM TRAP & CHK VLV)	3522, 3523, 3525 & 3533
CSP070	<i>No OEM Part # Available</i>	PM PKG (STM TRAP & CHK VLV)	3622 & 3633
CSP041	61301602510	SPARE PARTS PACKAGE	3525
CSP042	61301602498	SPARE PARTS PACKAGE	3522
CSP043	61301602499	SPARE PARTS PACKAGE	3523
CSP044	61301602500	SPARE PARTS PACKAGE	3533
CSP045	61301602501	SPARE PARTS PACKAGE	3522
CSP046	61301602502	SPARE PARTS PACKAGE	3523
CSP047	61301602503	SPARE PARTS PACKAGE	3533
NEW! CSP110	61301605861	PM PACK (LOWER)	500 Series
NEW! CSP111	61301605862	PM PACK (UPPER)	400 Series & 500 Series
NEW! CSP112	61301605863	PM PACK (LOWER)	400 Series
NEW! CSP113	61301605864	PM PACK (DOUBLE DOOR)	400 Series & 500 Series

Explore our complete line of parts to fit **Getinge/Castle Bulk Sterilizers**

including Door Gaskets, Solenoid Vales and Repair Kits, Steam Traps, Check Valves and Repair Kits, and more compatible with a variety for models.

Click or scan
to learn more

