



New Line of Parts to fit


Medivators®

DSD-201® & DSD Edge® AER Systems!


STERILIZERS - MINNTECH (MEDIVATORS)




RPI Part #MTA053
OEM Part # (No OEM Part # Available)
MANIFOLD ASSEMBLY
(AIR VALVE, 8 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge




RPI Part #MTM054
OEM Part # (No OEM Part # Available)
MANIFOLD (AIR VALVE, 8 STN)
 • Includes all parts as shown
 • Included in: Manifold Assembly (Air Valve, 8 STN) (RPI Part #MTA053)
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge




RPI Part #MTA034
OEM Part #MM01-0035
MANIFOLD ASSEMBLY
(AIR VALVE, 6 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge




RPI Part #MTM035
OEM Part # (No OEM Part # Available)
MANIFOLD (AIR VALVE, 6 STN)
 • Includes all parts as shown
 • Included in: Manifold Assembly (Air Valve, 6 STN) (RPI Part #MTA034)
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge




RPI Part #MTA036
OEM Part #MM01-0028
MANIFOLD ASSEMBLY
(AIR VALVE, 4 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge



RPI Part #MTM037
OEM Part # (No OEM Part # Available)
MANIFOLD (AIR VALVE, 4 STN)
 • Includes all parts as shown
 • Included in: Manifold Assembly (Air Valve, 4 STN) (RPI Part #MTA036)
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge



RPI Part #RXM047
OEM Part # (No OEM Part # Available)
PURGE MUFFLER (SINTERED)
 • Material: Sintered bronze and nickel plated steel
 • Port: 1/8" MP
 • Size: 7/16" Hex x 1-1/8" lg.
 • 40 micron element
Fits: Air Valve Manifold
Models: DSD-201 & DSD Edge



RPI Part #MTV030
OEM Part #41600-109
SOLENOID VALVE (MAC)
 • Includes all parts as shown
Fits: Air Delivery System
Models: DSD-201 & DSD Edge


Also available: Tubing (.106" ID x 5/32" OD Black) (RPI Part #RXT010) and Tubing (.106" ID x 5/32" OD Orange) (RPI Part #RXT011)




replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@riparts.com • WEBSITE www.riparts.com


The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Flyer Code 1606MT Rd2 Show 05/20



RPI Part #MTA038
OEM Part # (No OEM Part # Available)
MANIFOLD ASSEMBLY (8 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge




RPI Part #MTM039
OEM Part #MM01-0021
MANIFOLD (8 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis, Manifold Assembly
Models: DSD-201 & DSD Edge




RPI Part #MTM051
OEM Part #75227-516
SUB MANIFOLD (CHECK VALVES)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis, Manifold Assembly
Models: DSD-201 & DSD Edge




Also available: Collet Release Wrench Set (RPI Part #RXT004), Wrench (Low Profile) (RPI Part #RXT003), Tubing (.106" ID x 5/32" OD Black) (RPI Part #RXT010), Tubing (.106" ID x 5/32" OD Orange) (RPI Part #RXT011), and Locking Clips (RPI Part #RXC057)



RPI Part #MTA040
OEM Part #MM01-0030
MANIFOLD ASSEMBLY (6 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge



RPI Part #MTM041
OEM Part #MB02-0551
MANIFOLD (6 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis, Manifold Assembly
Models: DSD-201 & DSD Edge





RPI Part #MTM051
OEM Part #75227-516
SUB MANIFOLD (CHECK VALVES)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis, Manifold Assembly
Models: DSD-201 & DSD Edge




RPI Part #MTV027
OEM Part #78401-260


CHECK VALVE
 • Includes all parts as shown
Fits: Fluid Manifold
Models: DSD-201 & DSD Edge




RPI Part #MTA042
OEM Part #ASM1-0086
MANIFOLD ASSEMBLY (2 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge



RPI Part #MTM043
OEM Part #MB02-0552
MANIFOLD (2 STN)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis
Models: DSD-201 & DSD Edge



RPI Part #MTM052
OEM Part #MB02-0549
SUB MANIFOLD (CHECK VALVES)
 • Includes all parts as shown
 • See **Also available** below
Fits: Chassis, Manifold Assembly
Models: DSD-201 & DSD Edge



RPI Part #MTH026
OEM Part #43580-100
WIRING HARNESS
 • Length: 18" lg.
 • Also available: Solenoid Valve (Alcohol/Detergent) (RPI Part #MTV020)
 • Included in: Manifold Assembly (8 STN) (RPI Part #MTA038) and Manifold Assembly (6 STN) (RPI Part #MTA040)
Fits: Alcohol/Detergent Valves
Models: DSD-201 & DSD Edge



IMPORTANT NOTE: The OEM Part #41600-124 refers to two different valve block configurations. Valve Block Assembly (RPI Part #MTA016) includes (2) elbow fittings. Valve Block Assembly (RPI Part #MTA018) includes (1) elbow fitting and (1) straight fitting. Verify the configuration needed for your application.

ALSO AVAILABLE: Valve Block Rebuild Kit (RPI Part #MTK055), Tubing (.106" ID x 5/32" OD Black) (RPI Part #RXT010), Tubing (.106" ID x 5/32" OD Orange) (RPI Part #RXT011), and Wrench (Low Profile) (RPI Part #RXT003)

RPI Part #MTA016
OEM Part #41600-124




VALVE BLOCK ASSEMBLY

- Material: Hard anodized aluminum body
- Includes all parts as shown
- Included in: Manifold Assembly (2 STN) (RPI Part #MTA042)
- See **Important Notes** above

Fits: Manifold Assembly

Models: DSD-201 & DSD Edge



RPI Part #MTA017
OEM Part #41600-123




VALVE BLOCK ASSEMBLY

- Material: Hard anodized aluminum body
- Includes all parts as shown
- Included in: Manifold Assembly (RPI Part #'s MTA038 & MTA040)
- See **Important Note** above

Fits: Manifold Assembly

Models: DSD-201 & DSD Edge



RPI Part #MTA018
OEM Part #41600-124




VALVE BLOCK ASSEMBLY

- Material: Hard anodized aluminum body
- Includes all parts as shown
- Included in: Manifold Assembly (RPI Part #'s MTA038 & MTA040)
- See **Important Note** above

Fits: Manifold Assembly

Models: DSD-201 & DSD Edge



RPI Part #MTK055
OEM Part # (No OEM Part # Available)




VALVE BLOCK REBUILD KIT

- Includes all parts as shown
- Retaining ring pliers are necessary to remove and reinstall the internal retaining ring when rebuilding the Valve Block Assemblies (RPI Part #'s MTA016, MTA017 & MTA018).

Fits: Valve Block Assembly

Models: DSD-201 & DSD Edge



RPI Part #RXH083
OEM Part #MS05-0055




SCREW (#8-32 X 2-3/8" SOC HD)

- 8 per package; Material: Stainless steel
- The OEM valve block assembly's mounting screws are too short to use in RPI's manifolds, these longer screws must be used
- Included In: Valve Block Assembly (RPI Part #'s MTA016, MTA017 & MTA018)

Fits: Valve Block Assembly

Models: DSD-201 & DSD Edge



RPI Part #MTV020
OEM Part #MV01-034




VALVE ASSEMBLY (ALCOHOL/DETERGENT)

- Manifold mount; Manual operator
- 120VAC / 50-60 Hz., 8 Watt (240 Ω)
- Viton® seals
- Includes: Coil (RPI Part #MTC023), Solenoid Valve Repair Kit (RPI Part #MTK021), (1pc) Valve Seal (RPI Part #MTS022) and mounting hardware
- Also available: Wiring Harness (Solenoid Valve) (RPI Part #MTH026)
- Included in: Manifold Assembly (8 STN) (RPI Part #MTA038) and Manifold Assembly (6 STN) (RPI Part #MTA040)

Fits: Manifold Assembly

Models: DSD-201 & DSD Edge



RPI Part #MTC023
OEM Part #MV01-0032




COIL (SOLENOID VALVE)

- 120 VAC 50/60 Hz., 8W (240 Ω)

Fits: Alcohol/Detergent Valves

Models: DSD-201 & DSD Edge



RPI Part #MTK021
OEM Part # (No OEM Part # Available)




SOLENOID VALVE REPAIR KIT

- Includes: Plunger, spring, (1 pc) Seal (RPI Part #MTS022), O-Ring, bonnet and pal nut

Fits: Alcohol/Detergent Valves

Models: DSD-201 & DSD Edge



RPI Part #MTV050
OEM Part #MV01-0039




CHECK VALVE

- Cracking Pressure: 5 PSI
- Ports: 1/4" Push-In
- Material: White Acetal; Size: .59" dia. x 1.82" lg.

Fits: Alcohol and Detergent Lines

Models: DSD-201 & DSD Edge



RPI Part #MTV056
OEM Part #MV01-0027




CHECK VALVE

- Cracking Pressure: 10 PSI
- Ports: 1/4" MPT
- Material: Brass body, Viton® seal; Size: 9/16 Hex x 1.13" lg.

Fits: Air Tank

Models: DSD-201 & DSD Edge



RPI Part #RXV046
OEM Part #MV01-0038




PRESSURE RELIEF VALVE

- Cracking Pressure: 75 PSI
- Material: Brass; Port: 1/4" MPT

Fits: Air Tank

Models: DSD-201 & DSD Edge





ADDRESS SERVICES REQUESTED

New Line of Parts to fit Medivators® DSD-201® & DSD Edge® AER Systems!

RPI Part #RPF985
OEM Part #ME01-0033




ELBOW FITTING
(5/32" PUSH-IN x #10-32 MALE)

- 5 per package
- Material: Nickel plated brass
- Includes: Sealing O-Ring
- Included in: Valve Block Assembly (RPI Part #'s MTA016 and MTA018)

Fits: #10-32 Female Ports

Models: DSD-201 & DSD Edge



RPI Part #RPF986
OEM Part # (No OEM Part # Available)




MALE CONNECTOR
(5/32" PUSH-IN x #10-32 MALE)

- 5 per package
- Material: Nickel plated brass
- Includes: Sealing O-Ring
- Included in: Valve Block Assembly (RPI Part #'s MTA016 and MTA018)

Fits: #10-32 Female Ports

Models: DSD-201 & DSD Edge




RPI Part #RXF050
OEM Part #MA01-0076
FITTING (3/8" STEM x 1/4" MPT)



- 2 per package
- Material: White Polypropylene
- Included in: Manifold Assembly (8 STN) (RPI Part #MTA038), Manifold (8 STN) (RPI Part #MTA039), Manifold Assembly (6 STN) (RPI Part #MTA040), and Manifold (6 STN) (RPI Part #MTA041)
- See **Also available** to the right

Fits: 3/8" OD Tubing and 1/4" FPT Port

Models: DSD-201 & DSD Edge



RPI Part #RXF045
OEM Part #ME01-0026



ELBOW FITTING
(3/8" PUSH-IN x 3/8" STEM)

- 2 per package
- Material: White Polypropylene
- See **Also available** to the right

Fits: 3/8" OD Tubing

Models: DSD-201 & DSD Edge



RPI Part #RXF051
OEM Part #ME01-0027




ELBOW FITTING
(3/8" PUSH-IN x 3/8" PUSH-IN)

- 2 per package
- Material: White Polypropylene
- See **Also available** to the right

Fits: 3/8" OD Tubing

Models: DSD-201 & DSD Edge



RPI Part #RXF055
OEM Part # (No OEM Part # Available)




FITTING
(1/4" PUSH-IN X 3/8" STEM)

- 2 per package
- Material: White Polypropylene
- See **Also available** to the right

Fits: 1/4" OD Tubing

Models: DSD-201 & DSD Edge



RPI Part #RXF056
OEM Part #MT02-0017




TEE FITTING (3/8" PUSH-IN)

- Material: White Polypropylene
- See **Also available** to the right

Fits: 3/8" OD Tubing

Models: DSD-201 & DSD Edge



RPI Part #MTR049
OEM Part # (No OEM Part # Available)




RESTRICTOR

- Used to test drain times and simulate the pressure of an active scope.

Fits: Basin Fluid Inlet Line

Models: DSD-201 & DSD Edge



Also available ...

RPI Part #RXC057
OEM Part #ML05-0001



LOCKING CLIPS

- 10 per package
- Material: Red Acetal; Size: 3/8"

Fits: Used to Lock Collets on John Guest 3/8" Push-In Fittings

Models: DSD-201 & DSD Edge




RPI Part #RXT004
OEM Part #47049-157

COLLET RELEASE WRENCH SET

- John Guest collet wrench set
- Sizes: 3/16", 1/4", 5/16", 3/8" and 1/2"
- Also available: Locking Clips (RPI Part #RXC057)

Fits: John Guest Push-In Fittings

Models: DSD-201 & DSD Edge



RPI Part #RXT003
OEM Part #MW03-0002

WRENCH (LOW PROFILE)

- Material: Stainless steel
- Double Ended Flat Wrench: 7/16" x 10mm
- Used to secure the valve stem flat on Medivators pneumatic valve.
- Included in: PM Kit (RPI Part #MTK009) and PM Kit (RPI Part #MTK033)

Models: DSD-201 & DSD Edge

