

**New  
Parts  
to fit**

## Midmark® • Ritter® Exam Tables and Chairs!



**RPI Part #MIA315**  
**OEM Part #016-0358-00/**  
**002-0391-00/002-0564-10**



**ACTUATOR (BACK & TILT)**

- 2.2A @ 115VAC
- Includes: (2 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (2 pcs) Cable Tie (14" Black) (RPI Part #RPT816), and (3) QC terminal (1/4" male, fully insulated)
- Also available: Capacitor (36-43µf; 330VAC) (RPI Part #MIC219) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Seat Section

**Models:** 414, 416 & 417 (see Note below)



**RPI Part #MIA316**  
**OEM Part #016-0338-00/**  
**002-0296-00/002-0565-10**



**ACTUATOR (BASE)**

- 3.2A @ 115VAC
- Includes: (4 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (1 pc) Threadlocker 262 (RPI Part #RXA007), and (3) QC terminal (1/4" male, fully insulated)
- Also available: Capacitor (50µf, 370VAC) (RPI Part #MIC321) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Base Section

**Models:** 414 & 417 (see Note below)



**RPI Part #MIC321**  
**OEM Part #015-0438-06/**  
**002-1058-00**




**CAPACITOR (50µf, 370 VAC, 50/60 Hz)**

- Run capacitor
- Size: 2.00" dia. x 3.75" lg.
- Quad 1/4" Tabs
- Includes mounting bracket

**Fits:** Back Section

**Models:** 414 & 417 (see Note below)



**RPI Part #MIA317**  
**OEM Part #002-0727-00**



**ACTUATOR (BACK & TILT)**

- 2.2A @ 115VAC
- Includes: (2) Limit Switch (RPI Part #MIS175), (2 pcs) Cable Tie (8" Grey) (RPI Part #RPT397), (4 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (2 pcs) Wire Nut (RPI Part #RPT017), (3 pcs) QC Terminal (1/4" Female) (RPI Part #RPT472), and pre-installed modified wire harness
- Also available: Capacitor (36-43 µf; 330VAC) (RPI Part #MIC219) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Seat Section

**Models:** 416 & 417 (see Note below)



**RPI Part #MIA313**  
**OEM Part #016-0214-00/**  
**002-0596-02/002-0437-00**



**ACTUATOR (BASE)**

- 2.2A @ 115VAC
- Includes: (2 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (2 pcs) Cable Tie (14" Black) (RPI Part #RPT816), (3) #6 spade terminal, and (1 pc) Threadlocker 262 (RPI Part #RXA007)
- Also available: Capacitor (41-53µf; 330VAC) (RPI Part #MIC183) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Base Section

**Model:** 413 & 415 (see Note below)



**RPI Part #MIA314**  
**OEM Part #002-0596-02/**  
**016-0213-00/002-0435-00**



**ACTUATOR (SEAT)**

- 2.2A @ 115VAC
- Includes: (2 pcs) Cable Tie (14" Black) (RPI Part #RPT816), (2 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (3) QC terminal (1/4" male, fully insulated), (3) QC terminal (1/4" female, fully insulated), and (1 pc) Threadlocker 262 (RPI Part #RXA007)
- Also available: Capacitor (41-53µf; 330VAC) (RPI Part #MIC183), Limit Switch (RPI Part #MIS175) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Seat Section

**Models:** 413 & 415 (see Note below)

**Note:** Not all parts fit all Serial #'s for the models stated. Please see [www.rpiparts.com](http://www.rpiparts.com) for Serial # breakdown of the models these parts fit.



**RPI Part #MIA311**  
**OEM Part #016-0230-00/  
002-0344-00/029-0682-00**



**ACTUATOR (BASE)**

- 2.2A @ 115VAC
- Includes: (2 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (2 pcs) Cable Tie (14" Black) (RPI Part #RPT816), and (3) QC terminal (1/4" male, fully insulated)
- Also available: Capacitor (36-43µf; 330VAC) (RPI Part #MIC219) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Base Section  
**Model:** 405 (see Note below)




**RPI Part #MIA312**  
**OEM Part #016-0331-00/  
002-0345-00/029-0683-00**



**ACTUATOR (BACK)**

- 1.5A @ 115VAC
- Includes: (2 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (2 pcs) Cable Tie (14" Black) (RPI Part #RPT816), and (3) QC terminal (1/4" male, fully insulated)
- Also available: Capacitor (36-43µf; 330VAC) (RPI Part #MIC219) and Brake & Coupler Kit (RPI Part #MIK167)

**Fits:** Back Section  
**Model:** 405 (see Note below)



**RPI Part #MIA274**  
**OEM Part #002-0627-00**



**ACTUATOR (BASE)**

- 3.2A @ 115VAC
- Includes: Limit Switch (RPI Part #MIS175), (4 pcs) Cable Tie (12" Black) (RPI Part #RPT375), (1 pc) Cable Tie (8" Black) (RPI Part #RPT279), and (1 pc) Threadlocker 242 (RPI Part #RPA032)
- Also available: Capacitor (35µf; 370VAC, 50/60 Hz) (RPI Part #MIC220)

**Fits:** Base Section  
**Model:** 419 (see Note below)




**RPI Part #MIC318**  
**OEM Part #9A93001**



**HAND CONTROL**

- Includes: Hand Control Panel (RPI Part #MIP319), Hand Control PC Board (RPI Part #MIB320), Coiled Cord (RPI Part #MIC174), Top End Cap (RPI Part #MIC172), Bottom End Cap (RPI Part #MIC173), and pressure sensitive adhesive velcro strip

**Fits:** Chassis Connector  
**Models:** 414, 416 & 417 (see Note below)



**RPI Part #MIB320**  
**OEM Part #015-0568-00**



**HAND CONTROL PC BOARD**

- Included in: Hand Control (RPI Part #MIC318)
- Also available: Hand Control Panel (RPI Part #MIP319) and Coiled Cord (RPI Part #MIC174)

**Fits:** Hand Control  
**Models:** 414, 416 & 417 (see Note below)



**RPI Part #MIP319**  
**OEM Part #015-0546-00**

**HAND CONTROL PANEL**

- Size: 2.16" x 6.61"; Metal backing plate
- Included in: Hand Control (RPI Part #MIC318)
- Also available: Hand Control PC Board (RPI Part #MIB320), Coiled Cord (RPI Part #MIC174), Top End Cap (RPI Part #MIC172), and Bottom End Cap (RPI Part #MIC173)

**Fits:** Hand Control  
**Models:** 414, 416 & 417 (see Note below)



**RPI Part #MIS325**  
**OEM Part #015-0543-00/015-1587-00**



**SWITCH**

- Rocker switch (On-Off), Legend: I/O
- 16A @ 125VAC; 10A @ 250VAC; SPST
- 1/4" Tabs
- Panel cutout: 1.125" x .550", .025" - .187" Thk panel

**Fits:** Seat Section  
**Model:** 414, 416 & 417



**RPI Part #MIS326**  
**OEM Part #015-0635-00**



**LIMIT SWITCH**

- 11A @ 125/250 VAC; SPDT
- 3/16" Tabs

**Fits:** Base Section  
**Models:** 414 & 417 (see Note below)



**RPI Part #MIC294**  
**OEM Part #015-1000-00/**  
**015-1293-00/015-1396-00/**  
**015-1342-00**



**MAIN COILED HARNESS**

- Connects main PC board to back and tilt actuators.
- Material: 13 conductor cable; black jacket; all connectors and terminals pre-installed
- Includes: (4 pcs) loop clamps
- Depending on the model, not all connectors will be used
- Also available: Sensor Coiled Harness (RPI Part #MIC297)

**Fits:** Base Section

**Models:** 411, 419 & Evolution 75 (see Note below)

**RPI Part #RXM085**  
**OEM Part #015-1653-00/**  
**002-10454-00**




**POWER ENTRY MODULE**

- 10A @ Selectable 115 VAC/230 VAC, 3/16" tabs
- For use with C13 power cords
- Fuses not included, accepts up to 2 ea 1/4" dia. x 1-1/4" lg. or 5mm x 20mm fuses
- Panel cutout: 1.06" x 2.20", .031" - .079" panel Thk

**Fits:** Chassis

**Models:** 230, 630, 641, 647 & Evolution 75 (see Note below)

**RPI Part #MIC297**  
**OEM Part #015-1118-00**



**SENSOR COILED HARNESS**

- Connects main PC board to back and tilt position sensors.
- Material: 7 conductor cable; black jacket; all connectors pre-installed
- Includes: (2 pcs) loop clamps
- Also available: Main Coiled Harness (RPI Part #MIC294)

**Fits:** Base Section

**Models:** 411, 419 & Evolution 75 (see Note below)

**RPI Part #RXS113**  
**OEM Part #015-0007-04**



**HARNESS STRAP**

- 6 per package
- Size: 2-3/4" long, black nylon
- Mounts in a 1/4" dia. hole x 3/32" Thk max panel
- Use to strap cables, wires, and water and air supply lines

**Fits:** Chassis

**Models:** 222, 223, 622 & 623 (see Note below)

**RPI Part #MIC339**  
**OEM Part #015-0066-11/**  
**015-0066-04**



**HOSPITAL GRADE POWER CORD (15A @ 125VAC, 10 ft.)**

- 10 ft. lg. with jacket stripped 8"; Gray
- 12/3 SJT, 105° C; 15A @ 125 VAC
- Plug Type: NEMA 5-15 hospital grade
- Includes: (1 pc) Strain Relief Bushing (RPI Part #RXB112), (4 pcs) Cable Tie (8" Black) (RPI Part #RPT279), (6 pcs) Cable Tie (8" Black) (RPI Part #RPT085), connector housing with (3) crimp terminal, (1) #10 ring terminal, and (2) #6 spade terminal

**Fits:** Chassis

**Models:** 411, 414, 416 & 417

**RPI Part #RXB112**  
**OEM Part #015-0002-05**  
**Industry Part #SR 31-2**



**STRAIN RELIEF BUSHING**

- 3 per package
- Accommodates .425" - .475" dia. round cord
- Mounting hole: 7/8" dia.
- Included in: Hospital Grade Power Cord (15A @ 125VAC, 10 ft.) (RPI Part #MIC339)

**Fits:** Chassis

**Models:** 411, 416 & 417

**ALSO AVAILABLE:**

• **HOSPITAL GRADE POWER CORD (15A @ 125VAC, 10 ft.)**  
**RPI Part #RPC260**  
**OEM Part # (No OEM Part # Available)**

• **FUSE (6.3A, 250V)**  
**RPI Part #RXF036**  
**OEM Part #015-0346-40**



• **STRAIN RELIEF BUSHING PLIERS**  
**RPI Part #RPT172**  
**OEM Part # (No OEM Part # Available)**

• **HOSPITAL GRADE POWER CORD (15A @ 125VAC, 8 ft.)**  
**RPI Part #RPC554**  
**OEM Part # (No OEM Part # Available)**

• **FUSE HOLDER**  
**RPI Part #RPH674**  
**OEM Part #015-1259-02**

• **GAS CYLINDER**  
**RPI Part #MIC298**  
**OEM Part #016-0356-00**

# NEW ACTUATORS, HAND CONTROL, AND MORE!




## New Parts to fit Midmark® • Ritter® Exam Tables and Chairs!

ADDRESS SERVICES REQUESTED

**replacement parts industries, inc.**  
 625 Cochran Street • Simi Valley, California • 93065



RESORTED STANDARD  
 U.S. POSTAGE PAID  
 VAN NUYS, CA  
 PERMIT NO. 2399




**RPI Part #MIB322**  
**OEM Part #015-0568-01/002-0347-05**

### FOOT CONTROL PC BOARD

- Also available: Coiled Cord (RPI Part #MIC174)

**Fits:** Foot Control

**Models:** 414, 416 & 417



**RPI Part #MIK336**  
**OEM Part #002-0868-02**

### SLIDE KIT

- Replacement slides for the foot extension
- Includes: Two slides, all necessary mounting hardware and (1 pc) Threadlocker (242) (RPI Part #RPA032)

**Fits:** Seat Section Weldment

**Models:** 414, 416 & 417