

New line of parts to fit ...

STERIS® AMSCO® V-PRO®

Low Temperature Sterilization Systems

RPI Part #SSP058
OEM Part #P764333-354



PM PACK (GENERAL)

- Includes: HEPA Filter (RPI Part #SSF059), (4) Exhaust Filter (RPI Part #SSF060), (1 pc) Male Nut (1/8" OD Tubing) (RPI Part #RXF073), (1 pc) Ferrule (1/8" OD Tubing) (RPI Part #RXF074), Tee (1/8" OD Tubing) (RPI Part #RXF075), Battery (Lithium) (RPI Part #RXB076), Air Intake Filter (RPI Part #CMF282), Filter Element (RPI Part #SSE061), Solenoid Valve Repair Kit For 1/4" Port Valves (RPI Part #AMK175), and (2) Isolation Valve (3-Way) (RPI Part #SSV062)
 - Also available: PM Pack (Vacuum) (RPI Part #SSP065)
- Models:** V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSP065
OEM Part #P764334-481



PM PACK (VACUUM)

- Includes: Air Intake Filter (Metal Style) (RPI Part #CMF002), Exhaust Filter (RPI Part #SSF060), Filter Element (RPI Part #SSE066), (2) Cartridge Filter (RPI Part #SSF067), (1 pc) Disposal Bag (RPI Part #RPB838), (5 pcs) Absorbent Pad (RPI Part #RPP958), (3) Vacuum Pump Oil (Leybold Pump) (RPI Part #RXL034), and Drain Plug Gasket (RPI Part #SSG068)
 - All parts also sold separately
 - Also available: PM Pack (General) (RPI Part #SSP058)
- Models:** V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSP063
OEM Part #P764335-558



PM PACK (GENERAL)

- Includes: HEPA Filter (RPI Part #SSF059), (1 pc) Male Nut (1/8" OD Tubing) (RPI Part #RXF073), (1 pc) Ferrule (1/8" OD Tubing) (RPI Part #RXF074), and (1 pkg) Filter (RPI Part #SSF064).
 - All parts also sold separately
- Model:** V-PRO 60





RPI Part #SSV062
OEM Part #P093940-125

Isolation Valve (3-Way)

- 24VDC, 1.5W; 19.70" lg.
 - #24 AWG wires
 - Ports: 1/4-28; 3-Way
 - Included in: PM Pack (General) (RPI Part #SSP058)
- Fits:** Injection Assembly
Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSW072
OEM Part #P129603-014

SEALING WASHER

- 2 per package
 - Material: Red Silicone
 - Size: .290" ID x .688" OD x .093" thk
- Fits:** Oil Mist Eliminator
Models: V-PRO 60



x2

RPI Part #RXF075
OEM Part #P387348-974

TEE (1/8" OD TUBING)

- 1/4-28 for 1/8" OD tubing fitting
 - Material: PEEK®
 - Included in: PM Pack (General) (RPI Part #SSP058)
- Fits:** Hydrogen Peroxide Delivery System
Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



RPI Part #RXL034
OEM Part #P764334-644, P764335-529, P764335-529

VACUUM PUMP OIL (LEYBOLD PUMP)

- 1 liter (32 fl. oz.) bottle
 - Molykote L-1668 FM oil
 - Includes: Yorker spout
 - Included in: PM Pack (Vacuum) (RPI Part #SSP065)
 - Also available: Bottle (1 liter/32 fl. oz.) (RPI Part #RPB939), Disposal Bag (RPI Part #RPB838) and Absorbent Pad (RPI Part #RPP958)
- Fits:** Vacuum Pump
Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



RPI Part #RXF073
OEM Part #P387349-030

MALE NUT (1/8" OD TUBING)

- 6 per package
 - 1/4-28 x .61" lg for 1/8" OD tubing
 - Material: PEEK®
 - Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)
- Fits:** Hydrogen Peroxide Delivery System
Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



x6

RPI Part #AMK175
OEM Part #P910009-024

SOLENOID VALVE REPAIR KIT FOR 1/4" PORT VALVES


- For 1/4" port valves
 - Includes all parts as shown
 - Included in: PM Pack (General) (RPI Part #SSP058)
- Fits:** Solenoid Valve
Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #RXF074
OEM Part #P387349-031

FERRULE (1/8" OD TUBING)

- 6 per package
 - For 1/8" OD tubing
 - Material: Tefzel®
 - Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)
- Fits:** Hydrogen Peroxide Delivery System
Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



x6

RPI Part #SSD071
OEM Part #P129603-002

DIAPHRAGM (VALVE UMBRELLA)

- Material: Blue Silicone
 - Size: 18mm (.709") OD x 14.3mm (.563") lg.
- Fits:** Oil Mist Eliminator
Model: V-PRO 60



RPI Part #RXB076
OEM Part #P764332-350

BATTERY (LITHIUM)

- 3VDC Lithium
 - Size: .62" dia x 1.00" lg
 - Caution: Do not use if the battery is crushed or punctured.
 - Contact your state and local agencies for hazardous waste to determine proper disposal methods.
 - Included in: PM Pack (General) (RPI Part #SSP058)
- Fits:** Control Box Assembly
Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSG068
OEM Part #P764332-377

DRAIN PLUG GASKET

- Material: Viton®
 - Size: .615" ID x .875" OD x .080" thk
 - Included in: PM Pack (Vacuum) (RPI Part #SSP065)
- Fits:** Oil Mist Eliminator
Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #RXT077

OEM Part # (No OEM Part # Available)

TUBING (5/16" ID X 9/16" OD CLEAR)

- Sold by the foot
- Material: Industrial Grade PVC (Non-Food Grade)

Fits: Vacuum Pump Drain Use Only

Model: V-PRO 60



RPI Part #SSF067

OEM Part #P764334-369

CARTRIDGE FILTER

- 0.3 micron
- Size: 2.78" OD x 7.87" lg.
- Includes: Filter and (2) O-rings
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Fits: Oil Mist Eliminator

Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



RPI Part #SSF059

OEM Part #P101006-172

HEPA FILTER

- 1.2 micron
- Size: 2.50" OD x 7.60" lg.
- Ports: 3/8" FPT with 1/2" Barbs
- Included in: PM Pack (General) (RPI Part #'s SSP058 and SSP063)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



RPI Part #SSE061

OEM Part #P764333-215

FILTER ELEMENT

- 20 micron
- Size: .81" OD x .55" ID x 1.22" lg.
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #RPB838

OEM Part # (No OEM Part # Available)

DISPOSAL BAG

- 3 per package
- Used oil collection and disposal bag
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX




RPI Part #RPP958

OEM Part # (No OEM Part # Available)

ABSORBENT PAD

- 3 per package; Size: 9-1/4" x 14-1/2"
- Can be used for oil drip collection and spills
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



RPI Part #CMF002

OEM Part #P129385-485

AIR INTAKE FILTER (METAL STYLE)

- Size: 1.84" OD x 3.63" lg.
- Auto Muffler M02
- Included in: PM Pack (General) (RPI Part #SSP058)

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSF060

OEM Part #P129385-276

EXHAUST FILTER

- Size: .67" dia x 1.50" lg; 1/4" MPT
- Included in: PM Pack (General) (RPI Part #SSP058) and PM Pack (Vacuum) (RPI Part #SSP065)

Fits: Auto Drain Assembly

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSE066

OEM Part #P764334-368

FILTER ELEMENT

- 10 micron
- Size: 7.87" OD x 9.62" lg.
- Included in: PM Pack (Vacuum) (RPI Part #SSP065)

Fits: Air Intake

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #CMF282

OEM Part #P764333-087

AIR INTAKE FILTER

- 50 Microns
- Size: 2" dia. x 3.40" lg. x 1/4" MPT
- Includes: Air Filter Element (RPI Part #CME288)
- Included in: PM Pack (General) (RPI Part #SSP058)

Fits: Compressor

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #CME288

OEM Part # (No OEM Part # Available)

AIR FILTER ELEMENT

- 50 Microns; Size: 1.30" OD x .62" ID x 1.00" lg.
- Included in: Air Intake Filter (RPI Part #CMF282)

Fits: Air Intake Filter Housing

Models: V-PRO 1, V-PRO 1+ & V-PRO maX



RPI Part #SSF064

OEM Part #P129385-263

FILTER

- 2 per package
- Size: .30" dia x .90" lg; Port: M5 male
- Included in: PM Pack (General) (RPI Part #SSP063)

Fits: Injection Manifold (60), 3-way valve

Models: V-PRO 1, V-PRO 1+, V-PRO 60 & V-PRO maX



New PM Packs to fit STERIS® AMSCO® V-PRO® 1, V-PRO® 1+ & V-PRO® maX!

NEW



PM PACK (GENERAL)

RPI Part #SSP058

NEW



PM PACK (VACUUM)

RPI Part #SSP065

New Parts to fit STERIS® AMSCO® V-PRO® Sterilization Systems!

ADDRESS SERVICES REQUESTED

replacement parts industries, inc.
 625 Cochran Street • Simi Valley, California • 93065



PERMIT NO. 2399
 VAN NUYS, CA
 U.S. POSTAGE PAID
 PRESORTED STANDARD

New PM Pack to fit STERIS® AMSCO® V-PRO® 60!

NEW



PM PACK (GENERAL)

RPI Part #SSP063

Look inside for even more new parts to fit STERIS AMSCO V-PRO Low Temperature Sterilization Systems including:

- Battery (Lithium)
- Cartridge Filter
- Diaphragm (Valve Umbrella)
- Drain Plug Gasket
- Isolation Valve (3-Way)
- Exhaust Filter
- Filter Element
- Hepa Filter
- Sealing Washer
- Many More!