



RPI Part #ADK245
OEM Part #90.0593.00

FOOT CONTROL III SVC KIT


• Includes all parts as shown

Fits: Foot Control III (OEM Part #'s 38.1601.00, 38.1603.00, & 38.1700.00), 300 Series Foot Control (OEM Part #38.1821.00), & 500 Series Foot Control (OEM Part #38.1805.00)

Models: Components-Foot Control- Foot Control III, Foot Control 300 Series, & Foot Control 500 Series




RPI Part #ADC254
OEM Part #22.0050.00
SPRING CAP (FOOT CONTROL)



RPI Part #ADD022
OEM Part #38-0054-00
DIAPHRAGM




RPI Part #ADP083
OEM Part #22-0060-00
POPPET



RPI Part #ADP243
OEM Part #38.0760.00
PISTON (FOOT CONTROL III)




RPI Part #ADS081
OEM Part #10-0440-00
HELICAL COMPRESSION SPRING



RPI Part #ADS086
OEM Part #22-0580-00
HELICAL COMPRESSION SPRING




RPI Part #ADS087
OEM Part #22-0778-00
RELAY VALVE STEM



RPI Part #ADS091
OEM Part #013-011-00
HELICAL COMPRESSION SPRING



RPI Part #RPL090
OEM Part #98.0090.00
HIGH TEMP LUBRICANT



Shown Actual Size

RPI Part #RPO303
OEM Part #030.012.01
O-RING