



RPI Part #ADK257

OEM Part #38.1776.00 (Control Block Assy)

CONTROL BLOCK KIT (MINOR)


• Includes all parts as shown

Fits: Control Block 500 series

Models: Components-Control Block-500 series

(4x)


RPI Part #ADC236
OEM Part #38.1780.00



WATER CARTRIDGE ASSEMBLY

(2x)


RPI Part #ADC237
OEM Part #38.1783.00



WATER CARTRIDGE ASSY DRY

(1 PACK OF 5)


RPI Part #ADD238
OEM Part #38.1787.01



DIAPHRAGM

(2x)


RPI Part #ADR239
OEM Part #38.1800.00



BARB RETAINER

(1 PACK OF 5)


RPI Part #ADT242
OEM Part #38.1808.01 & 98.0867.02



STEM ADJUSTMENT TOOL

(4x)


RPI Part #ADS233
OEM Part #38.0516.00



STEM (WATER FLOW ADJ.)

(2x)

RPI Part #ADS234
OEM Part #38.0526.00



STEM (AIR COOLANT)

(4x)


RPI Part #ADS235
OEM Part #38.1785.01



STEM (DRIVE AIR)

RPI Part #RPC753
OEM Part # (No OEM Part # Available)


STORAGE CASE



(1 PACK OF 25)

RPI Part #RPH740
OEM Part #001.240.00

SOCKET HEAD CAP SCREW



(1 PACK OF 5)

RPI Part #RPL090
OEM Part #98.0090.00

HIGH TEMP LUBRICANT

(1 PACK OF 12)

RPI Part #RPO307
OEM Part #030.006.02

O-RING



Shown Actual Size

(1 PACK OF 12)

RPI Part #RPO484
OEM Part #030.003.01

O-RING



Shown Actual Size

(1 PACK OF 12)

RPI Part #RPO486
OEM Part #034.002.01

O-RING



Shown Actual Size

(1 PACK OF 12)

RPI Part #RPO752
OEM Part #035.034.01

O-RING



Shown Actual Size