





RPI Part #CSP041
OEM Part #61301602510

SPARE PARTS PACKAGE
• Includes all parts as shown
Fits: Machine in general
Models: 3525


RPI Part #CSP041 Spare Parts Package to fit Castle/Getinge® Bulk Sterilizers




RPI Part #AMK174
OEM Part #56161
SOLENOID VALVE REPAIR KIT




(2x)
RPI Part #AMK177
OEM Part #61301601739
SOLENOID VALVE REPAIR KIT



RPI Part #AMK183
OEM Part #54190
SOLENOID VALVE REPAIR KIT




RPI Part #CSD051
OEM Part #517153
DIAPHRAGM REPLACEMENT 1.0




RPI Part #CSF052
OEM Part #61301600303
FILTER, AIR-IN




RPI Part #CSG002
OEM Part #61301605017
DOOR GASKET (16" x 16")




(3x)
RPI Part #CSK048
OEM Part #68164
SOLENOID VALVE REPAIR KIT




RPI Part #CSK050
OEM Part #54187
SOLENOID VALVE REPAIR KIT



(3x)
RPI Part #CSK066
OEM Part #513369
CHECK VALVE REBUILD KIT




RPI Part #CSR028
OEM Part #86407
SOLID STATE RELAY




RPI Part #CSR029
OEM Part #86408
SOLID STATE RELAY




RPI Part #CSS054
OEM Part #45690
SCREEN, DETERGENT



RPI Part #CSV056
OEM Part #14010
CHECK VALVE - 1/4"



RPI Part #CSV057
OEM Part #19114
CHECK VALVE - 1/2"



RPI Part #CSV058
OEM Part #08045
CHECK VALVE - 3/4"

**Don't forget ...
RPI has other
PM and Spare Parts
Packages to fit
Castle/Getinge
Bulk Sterilizers.
[Click here for more.](#)**